



# Handle display

① Optional 12\*24mm/12\*36mm light spot size.

Different spot size meet different body parts hair removal needs .



spot size : 12\*24mm



spot size : 12\*36mm

② 10bars 1200w / 15bars 1800w. 10million—40million shots.  
Effective for permanent hair removal

③ 755nm and 808nm and 1064nm , suits for all skin and hair types.

# Cooling effect

The handle has a cooling effect that improve treatment comfort.



spot size : 12\*24mm



spot size : 12\*36mm


Waves 755nm 808nm 1064nm

755 ( / ) ( 808 ) ( 1064 )


1064nm 1064nm (1064nm 1064nm) 1064nm 1064nm  
1064nm 1064nm 1064nm 1064nm 1064nm

# working principle


- 755 nm for white skin (fine, golden hair removal );
- 808nm for yellow neutral skin (brown hair removal);
- 1064nm for black skin (black hair removal)



A natural




B energy absorbed by hair



C hair follicle died

# working principle

Triple wavelength laser spots work together



□□□□□□ □□□□□□ □□□□ □□□□

## Scope of application

Suitable for any unwanted hair on areas like face, arms, armpits, chest, back, bikini, legs...





# SHIPPING

Welcome to visit our factory



## In terms of following trade details:

- 1)The shipping cost are based on DHL cost to your country
- 2)Delivery time:3-5 days after receiving payments.
- 3)Payment terms:T/T(bank transfer),Western Union.
- 4)Warranty:2 years for the whole machine
- 5)Warranty Policy:one month guarantee replacement (we bear all shipping cost),2years warranty,lifelong maintenance.
- 6)Packing method: Aluminum box with hard plastic foam.
- 7)Free-training to customers:Free We could offer DVD and User's manual and 10 Hours online service to our customers
- 8)Certifications:CE TUV ISO 13485
- 9)MOQ:1 UNIT
- 10)The quotation here above is valid within 1 month