



EXFU·FIND S
EMSCULPT

25% more muscle,
30% less fat,
now all the way to your flanks

* 此設備係根據美國FDA批准之非侵入式電磁刺激技術，透過電磁波刺激腹部肌肉，達到鍛煉肌肉、減少脂肪之效果。此設備係根據美國FDA批准之非侵入式電磁刺激技術，透過電磁波刺激腹部肌肉，達到鍛煉肌肉、減少脂肪之效果。

* 此設備係根據美國FDA批准之非侵入式電磁刺激技術，透過電磁波刺激腹部肌肉，達到鍛煉肌肉、減少脂肪之效果。此設備係根據美國FDA批准之非侵入式電磁刺激技術，透過電磁波刺激腹部肌肉，達到鍛煉肌肉、減少脂肪之效果。



EASY TO OPERATE
**RANDOMLY ADJUST
MORE INTELLIGENT**
NEW APPEARANCE

* 此設備係根據美國FDA批准之非侵入式電磁刺激技術，透過電磁波刺激腹部肌肉，達到鍛煉肌肉、減少脂肪之效果。此設備係根據美國FDA批准之非侵入式電磁刺激技術，透過電磁波刺激腹部肌肉，達到鍛煉肌肉、減少脂肪之效果。

REASONABLE DESIGN

PRODUCT DETAILS

SAFE AND CONVENIENT



HANDLE

The handles support independent and synchronous work

POWER LOCK

Prevent the power supply from falling out, causing damage to the machine



SHELF

Put the tools you can use

STRAP STORAGE

The handles support independent and synchronous work



EASY TO OPERATE

HD 15.6 INCH LARGE SCREEN

NEW APPEARANCE

SCREEN DISPLAY

EASY TO OPERATE

ONE-BUTTON START

MORE CONVENIENT

COOLING SYSTEM

GOOD BATTERY LIFE



OPERATING PARTS


The location of the handle operation must be consistent with the interface selection location



INCREASE MUSCLE MASS

HIGH EFFICIENCY

UP TO 50.000 PULSES IN 30 MINUTE SESSION



VS



+16%

Estimated increase in
muscle mass

-19%

Estimated fat
reduction

4 CM (1.5 INCHES) LOSS IN WAISTLINE CIRCUMFERENCE

11% REDUCTION IN ABDOMINAL MUSCLE SEPARATION (CALLED DIASTASIS RECTI)

5 X INCREASE IN FAT METABOLISM (APOPTOTIC INDEX INCREASED FROM
19% TO 92% AFTER TREATMENT) 80% PATIENTS WITH VISIBLE LIFTING

EFFECT OF BUTTOCKS 96% PATIENT SATISFACTION


PRINCIPLE

HOW IT WORKS

EMSCULPT with HIFEM technology induces supramaximal muscle contractions not achievable by voluntary action. When exposed to supramaximal contractions, the muscle tissue is forced to adapt to such extreme condition.

The muscle responds with a deep remodeling of its inner structure, i.e., the growth of myofibrils (muscle hypertrophy) and creation of new protein strands and muscle fibers (muscle hyperplasia).

The process results in increased muscle density and volume.



Direct muscle stimulation, more efficient

More than 20,000 Sit-ups



Faster frequency

1-150Hz adjustable



More modes

5+1 optional



Higher efficiency

30 minutes treatment



Pays off faster

Speed up fat burning

DEVICE NAME:	EXFU TESLASCULPT MACHINE
MODEL NUMBER:	HI-EMT
APPLICATION:	FOR COMMERCIAL
SCREEN SIZE:	14INCH TRUE COLOR TOUCH SCREEN
PULSE DURATION:	300PS
FREQUENCY:	1-10HZ
SPOT SIZE:	2-10MM ADJUSTABLE
ENERGY:	20,000 SIT-UPS/50,000 MUSCLES CONTRACTIONS
AFTER-SALES SERVICE PROVIDED:	FREE SPARE PARTS, ONLINE SUPPORT, VIDEO TECHNICAL SUPPORT
WEIGHT:	65KGS
COOLING:	WATER COOLING+AIR COOLING
VOLTAGE:	AC220V/110V 50HZ/60HZ
DEMENSION OF MACHINE:	27*49*57CM
PACKAGE:	ALUMINUM BOX
WARRANTY:	2 YEARS

AMAZING EFFECT

Before & After



Before

Five weeks later

Eight weeks later



Before

After



Before

After



Before

After



Before

After