



 $^{^{*}}$

REASONABLE DESIGN

PRODUCT DETAILS

SAFE AND CONVENIENT



HANDLE

The handles support independent and synchronous work

POWER LOCK

Prevent the power supply from falling out, causing damage to the machine





SHELF

Put the tools you can use

STRAP STORAGE

The handles support independent and synchronous work



EASY TO OPERATE

HD 15.6 INCH LARGE SCREEN

NEW APPEARANCE

SCREEN DISPLAY

EASY TO OPERATE

ONE-BUTTON START

MORE CONVENIENT

COOLING SYSTEM

GOOD BATTERY LIFE

OPERATING PARTS


The location of the handle operation must be consistent with the interface selection location



INCREASE MUSCLE MASS


HIGH EFFCIENCY

UP TO 50.000 PULSES IN 30 MINUTE SESSION










PRINCIPLE

HOW IT WORKS

EMSCULPT with HIFEM technology induces supramaximal muscle contractions not achievable by voluntary actionWhen exposed to supramaximal contractions, the muscle tissue is forced to adapt to such extreme condition.

The musle responds with a deep remodeling of its inner structure, i.e, the growth of myofibr ils (muscle hypertrophy)and creation of new protein strands and muscle fibers (muscle hyperplasia).

he orocess resllpcreasecmoUscdensity andvolume





DEVICE NAME:	EXFU TESLASCULPT MACHINE
MODEL NUMBER:	HI-EMT
APPLICATION:	FOR COMMERCIAL
SCREEN SIZE:	14INCH TRUE COLOR TOUCH SCREEN
PULSE DURATION:	300PS
FREQUENCY:	1-10HZ
SPOT SIZE:	2-10MM ADJUSTABLE
ENERGY:	20,000 SIT-UPS/50,000 MUSCLES CONTRACTIONS
AFTER-SALES SERVICE PROVIDED:	FREE SPARE PARTS, ONLINE SUPPORT, VIDEO TECHNICAL SUPPORT
WEIGHT:	65KGS
COOLING:	WATER COOLING+AIR COOLING
VOLTAGE:	AC220V/110V 50HZ/60HZ
DEMENSION OF MACHINE:	27*49*57CM
PACKAGE:	ALUMINUM BOX
WARRANTY:	2 YEARS

AMAZING EFFECT

Before & After



Before

Five weeks later

Eight weeks later



Before

After

Before

After



Before



After



Before



After