

DIODE LASER



ICELLO



MLKJ es la mejor máquina de depilación láser comercial.

Según el principio de fototermodinámica selectiva, la longitud de onda, la energía y el ancho del pulso del láser se pueden ajustar adecuadamente para que el láser pase desde la superficie de la piel hasta la raíz del folículo piloso. La energía de la luz se absorbe y se convierte en energía térmica que destruye el tejido del folículo piloso, haciendo que el cabello no pueda regenerarse, lo que hace que los láseres MLKJ sean indoloros y permanentes con tecnología patentada.

La máquina de depilación láser es un método de tratamiento no invasivo, después del tratamiento no causará ningún daño especial a la piel circundante, no afectará la vida normal ni el trabajo.

□ Dispositivo certificado por la FDA y el CE médico

□ 3 en 1: tres tecnologías láser en un solo dispositivo

□ Láseres de Alejandrita, Diodo y Nd:YAG: tratan todo tipo de piel

- Fácil operación
- La capacitación en línea está disponible
- ¡50 millones de disparos por pieza de mano!

Pantalla de manija



1. Sin dolor. La punta de la pieza de mano se enfría hasta -32 grados Celsius y hace que el tratamiento sea cómodo y eficaz.
2. Resultados primarios del tratamiento. 2-3 sesiones de tratamiento pueden mostrar resultados permanentes
3. Prácticamente indoloro en contacto con el zafiro, comodidad y seguridad incomparables.
4. La punta de la cara aumenta la accesibilidad. Trate fácilmente áreas que normalmente son difíciles de alcanzar, incluidas las orejas, las fosas nasales y la glabella.



El precio de la máquina de depilación láser 755/808/1064nm.

>>808 single wave or 755nm/808nm/1064nm perfect triple wavelengths

>>America Imported Emitter, strong energy

>>10/12 pieces or 16 pieces Germany Dilas original laser bars, 1200W/1800W power is capable to meet all demands of hair removal on the whole body


Medical diode laser

50000000+shots
can treat 10000+clients

- *Not easy to be burned
- *Low maintenance cost
- *Stable quality



- >>1-10Hz fast frequency,faster, smooth and effective treatment
- >>Super sapphire touch cooling system ensures painless and comfortable treatment
- >> Italy imported stainless water pump, it can keep machine working long time without stop. Make sure the heat were removed very quickly
- >>Double filters to keep water quality



The real experience is visible

Through live test verification ,safe ,long,painless
and effective



Pantalla 4K de 14 pulgadas, experiencia táctil inteligente

- * 6 estilos de piel para elegir
- * Fácil de proporcionar habilidades operativas/servicios de capacitación.
- * El logotipo y el idioma gratuitos se pueden personalizar en la pantalla.



You'll get

Two major service guarantees
Worry-free after-sales


01 Professional after-sales group

24 hours offers the solution
Solve problems quickly with your active cooperation

02 Free maintenance

Warranty period, Send your replacement
parts for free

Fabricante de herramientas de belleza desde hace 15 años con certificación FDA, CE, ISO13485.
¡Estamos comprometidos a hacer que su negocio de belleza sea más fácil!



Establishment Registration & Device Listing

FDA Home Medical Devices Databases



New Search Back To Search Results

Establishment:
 WEIFANG MINGLIANG ELECTRONICS CO.,LTD.
 Floor 3, Building 4,Second Accelerator, Yuqing Eastern Street, High-Tech District
 Weifang Shandong, CN 261000
Registration Number: 3007549439
FEI Number*: 3007549439
Status: Active
Date Of Registration Status: 2022

Owner/Operator:
[Weifang Mingliang Electronics CO.,LTD.](#)
 Floor 3, Building 4,Second Accelerator, Yuqing Eastern Street, High-Tech District
 Weifang, Shandong CN 261000
Owner/Operator Number: 10082104

Official Correspondent:
 Pingping Lv
 Floor 3, Building 4,Second Accelerator, Yuqing Eastern Street, High-Tech District
 Weifang, Shandong CN 261000
Phone: 86-10-56335780

US Agent:
 Wen Xu
 4611 Maydew Mews Way Ellicott City MD
 Ellicott, MD US 21042
Phone: 443 8968583 Ext
Email: Infor@Believe-Med.Com

* Firm Establishment Identifier (FEI) should be used for identification of entities within the imports message set

Page Last Updated: 05/09/2022

Note: If you need help accessing information in different file formats, see Instructions for Downloading Viewers and Players.

Language Assistance Available: Español | 繁體中文 | Tiếng Việt | 한국어 | Tagalog | Русский | العربية | Kreyòl Ayisyen | Français | Polski | Português | Italiano | Deutsch | 日本語 | فارسی | English



Authentication

FDA,NQA13485,CB certification
and more Complete certification
materials



60000+ Customers

We have professional after sale
team and we provide 7*24
services



2 Years Warranty

2 years warranty for each
machine, All questions will get
reponse within 24 hours

¡Bienvenidos a la consulta! ¡Servicio en línea 7x24 horas!