

# DIODE LASER



# ICELLO



## **MLKJ est la meilleure machine d'épilation au laser commerciale.**

Selon le principe de la photothermodynamique sélective, la longueur d'onde, l'énergie et la largeur d'impulsion du laser peuvent être correctement ajustées afin que le laser passe de la surface de la peau au follicule pileux. L'énergie lumineuse est absorbée et convertie en énergie thermique qui détruit le tissu du follicule pileux, rendant les cheveux incapables de se régénérer, rendant les lasers MLKJ indolores et permanents grâce à une technologie brevetée.

La machine d'épilation au laser est une méthode de traitement non invasive. Après le traitement, elle ne causera aucun dommage particulier à la peau environnante et n'affectera pas la vie et le travail normaux.

□ Appareil certifié FDA et médical CE

□ 3 en 1 - trois technologies laser dans un seul appareil

□ Lasers alexandrite, diode et Nd:YAG - traitent tous les types de peau

- Opération facile
- Une formation en ligne est disponible
- 50 millions de clichés par pièce à main !

### **Affichage de la poignée**



1. Sans douleur. La pointe de la pièce à main refroidit jusqu'à -32 degrés Celsius et rend le traitement confortable et efficace.
2. Résultats primaires du traitement. 2-3 séances de traitement peuvent montrer des résultats permanents
3. Pratiquement indolore au contact du saphir, confort et sécurité inégalés
4. La pointe du visage augmente l'accessibilité. Traitez facilement les zones habituellement difficiles à atteindre, notamment les oreilles, les narines et la glabelle.



Le prix de la machine d'épilation au laser 755/808/1064nm.

>>808 single wave or 755nm/808nm/1064nm perfect triple wavelengths

>>America Imported Emitter, strong energy

>>10/12 pieces or 16 pieces Germany Dilas original laser bars, 1200W/1800W power is capable to meet all demands of hair removal on the whole body

# Medical diode laser

**50000000+shots**  
**can treat 10000+clients**

- \*Not easy to be burned
- \*Low maintenance cost
- \*Stable quality



- >>1-10Hz fast frequency,faster, smooth and effective treatment
- >>Super sapphire touch cooling system ensures painless and comfortable treatment
- >> Italy imported stainless water pump, it can keep machine working long time without stop. Make sure the heat were removed very quickly
- >>Double filters to keep water quality



# The real experience is visible

Through live test verification ,safe ,long,painless  
and effective



Écran 4K de 14 pouces, expérience tactile intelligente

- \* 6 styles de peau au choix
- \* Facile à fournir des compétences opérationnelles/des services de formation
- \* Le logo gratuit et la langue peuvent être personnalisés sur l'écran



# You'll get

Two major service guarantees  
Worry-free after-sales


## 01 Professional after-sales group

24 hours offers the solution  
Solve problems quickly with your active cooperation

## 02 Free maintenance

Warranty period, Send your replacement  
parts for free

Fabricant d'outils de beauté depuis 15 ans avec certification FDA, CE, ISO13485.  
Nous nous engageons à faciliter votre activité beauté !



# Establishment Registration & Device Listing

FDA Home Medical Devices Databases



New Search	Back To Search Results
<p><b>Establishment:</b> WEIFANG MINGLIANG ELECTRONICS CO.,LTD. Floor 3, Building 4,Second Accelerator, Yuqing Eastern Street, High-Tech District Weifang Shandong, CN 261000 <b>Registration Number:</b> 3007549439 <b>FEI Number*:</b> 3007549439 <b>Status:</b> Active <b>Date Of Registration Status:</b> 2022</p> <p><b>Owner/Operator:</b> <a href="#">Weifang Mingliang Electronics CO.,LTD.</a> Floor 3, Building 4,Second Accelerator, Yuqing Eastern Street, High-Tech District Weifang, Shandong CN 261000 <b>Owner/Operator Number:</b> <a href="#">10082104</a></p> <p><b>Official Correspondent:</b> Pingping Lv Floor 3, Building 4,Second Accelerator, Yuqing Eastern Street, High-Tech District Weifang, Shandong CN 261000 <b>Phone:</b> 86-10-56335780</p> <p><b>US Agent:</b> Wen Xu 4611 Maydew Mews Way Ellicott City MD Ellicott, MD US 21042 <b>Phone:</b> 443 8968583 Ext <b>Email:</b> <a href="mailto:Infor@Believe-Med.Com">Infor@Believe-Med.Com</a></p>	
<p>* Firm Establishment Identifier (FEI) should be used for identification of entities within the imports message set</p>	

Page Last Updated: 05/09/2022

Note: If you need help accessing information in different file formats, see Instructions for Downloading Viewers and Players.

Language Assistance Available: Español | 繁體中文 | Tiếng Việt | 한국어 | Tagalog | Русский | العربية | Kreyòl Ayisyen | Français | Polski | Português | Italiano | Deutsch | 日本語 | فارسی | English



### Authentication

FDA,NQA13485,CB certification  
and more Complete certification  
materials



### 60000+ Customers

We have professional after sale  
team and we provide 7\*24  
services



### 2 Years Warranty

2 years warranty for each  
machine, All questions will get  
reponse within 24 hours

**Bienvenue à la consultation ! Service en ligne 7x24 heures !**