



40K CAVITATION 5IN1

RF+VACUUM+CAVITATION
+LASER PADS SLIMMING MACHING

START



Advantage^e



Advantage

- **portable**

5 in 1 salon beauty system
saving shipping cost

5 handles

for slimming, skin tightening, wrinkle removal

- **Advantage**

high power cavitation 300W more effective

vacuum head help a lot for supus fat cell

2 years warranty

multi-color for option





Cavitation

Frequency:4Mhz

Power:150W

Burning explode fat cell



Monopolar RF

Frequency:40khz

Power:300W

Body shape fat reduction



Big tripolar RF

Frequency:4Mhz

Power:150W

Skin tightening lifting wrinkle removal



Small Tripolar RF

Frequency:1Mhz

Power:150W

Skin tightening lifting wrinkle removal



Vacuum

Frequency:1MHz

Pressure:-90kpa

Slimming Remove surplus fat cell

article	Machine amincissante cavitation vide RF
Écran LCD TFT	Écran tactile couleur de 8 pouces
Technologie	CavitationVideRF (RF monopolaireRF tripolaire) Tête de machine à cavitation RF à ultrasons Tête de liposuction sous vide
Pièces à main	Tête RF monopolaire Tête RF tripolaire (grande) Tête RF tripolaire (petite)
Fréquence de cavitation	40 kHz

Puissance maximale de cavitation	300w
RF tripolaire	1-10 MHZ (fréquence de sortie) 150 watts (puissance de sortie maximale)
RF monopolaire	1-10 MHZ (fréquence de sortie) 150 watts (puissance de sortie maximale)
Pression du vide	100 KPA
Taille	Taille de la machine : 47*31*38 cm Taille d'emballage du boîtier en aluminium: 54*43*46cm
Tension	220 V CA/110 V CA
N.W.	15KG
G. W.	25KG

Treatment before & after compare




Effet de l'application et du traitement :

- > Suitable for 6 different color skin, including the dark skin.
- > The sapphire temperature cools down 0~3°C, clients feel comfortable and painless in whole treatment.
- > Long continuous working time, ensure the works in Clinic/salon use.

AMAZING EFFECT

Before & After



Before


Five weeks later

Eight weeks later



Before

After



Before

After



Before

After



Before

After

Free logo service



[Consultez immédiatement!](#)