




Macro & Micro High Intensity Focused Ultrasound System
Non-Surgical Treatment Of Face Lifting, Body Firming And Body Contouring System!

Products Description




Use high energy ultrasound focusing on the treatment position, so that the skin tissue create the thermal, and make high-speed friction cells to stimulate collagen. Such heat effect will not hurt the epidermis, because the treatment fast and direct access to the treatment site within 0-0.5 seconds, without touching the surrounding tissue and can be passed directly to the superficial tendon membrane system (SMAS), so it can be done firm skin while pulling the muscle layer, thin face up to the progressive effect.




MF2.0MM

Unique ultra-thin
probe Clexis patented
technology

"Special cartridge for eyes treatment"



Body Treatment



• *Treatment Effect* •



| | |
|---------------------|---|
| Name | High Intensity Focused Ultrasound |
| Handles | Facial: 1.5mm, 3.0mm, 4.5mm |
| | 7D HF: MF2/6/9/13 L4-4.5 L7-1.5 L7-3.0 |
| Energy output | 0.1J-2.5J (0.1J per step) |
| Space between Point | 1.0mm-4.0mm |
| Shot length | 5mm-25mm |
| shots | Each cartridge 300000 shots |
| Screen Size | 15 inch |
| Voltage | 110V-240V, 50/60Hz |

Welcome to consult more product details!