





Macro & Micro High Intensity Focused Ultrasound System
 Non-Surgical Treatment Of Face Lifting, Body Firming And Body Contouring System!

Products Description




Use high energy ultrasound focusing on the treatment position, so that the skin tissue create the thermal, and make high-speed friction cells to stimulate collagen. Such heat effect will not hurt the epidermis, because the treatment fast and direct access to the treatment site within 0-0.5 seconds, without touching the surrounding tissue and can be passed directly to the superficial tendon membrane system (SMAS), so it can be done firm skin while pulling the muscle layer, thin face up to the progressive effect.




MF2.0MM

Unique ultra-thin
probe Clexis patented
technology

"Special cartridge for eyes treatment"



Body Treatment



• *Treatment Effect* •



Name	High Intensity Focused Ultrasound
Handles	Facial: 1.5mm, 3.0mm, 4.5mm
	7D HF: MF2/6/9/13
	L4-4.5 L7-1.5 L7-3.0
Energy output	0.1J-2.5J (0.1J per step)
Space between Point	1.0mm-4.0mm
Shot length	5mm-25mm
shots	Each cartridge 300000 shots
Screen Size	15 inch
Voltage	110V-240V, 50/60Hz

Welcome to consult more product details!