



Multifunctional elight ipl diode laser 808nm hair removal nd yag tatoo removal with rf

Model	Elight+808nm diode laser+nd yag Laser+RF
Application	1)Hair removal, freckle removal 2)Skin rejuvenation 3)Pigmentation treatment 4)Eliminate wrinkles 5)Acne clearance 6)Eyeline removal 7)Eyebrow embroidery removal 8)Tattoo removal 9)Freckles removal

Light Source	Intense Pulsed Light
Handpiece	One Elight handle One 808nm diode laser handle One nd yag laser handle and one RF handle
Wavelength	IPL:420/510/690~1200NM (IPL Standard) 585/640/755~1200nm(optional) SHR:640-690NM
Pulse Delay Time	1~50ms
Fluency	2~15ms
Screen	8.4 inch color touch screen
Spot size	12x20mm, 12x35mm, 15x50mm
Q-Switched Nd	ND: YAG laser
Wavelength	1064&532&1320nm
Pulse energy	0-800mj 0-300mj
Energy of SHR	1-50J
Cooling system	Aire and colse water circulation, treatment head with semiconductor cooling temperature adjustable
Power	1600-3000W
Electric source	-4-----0
Packing method	Aluminum Packing
Packing method	Aluminum Packing

Multifunctional elight ipl diode laser 808nm hair removal nd yag tatoo removal

808nm diode laser handle for super fast hair removal

e-light(IPL+RF) handle for skin rejuvenation ,skin tightening ,ance removal ,hair removal and skin care etc.

nd yag laser handle for tattoo removal ,eyeline removal etc

rf head for facial and body slimming ,skin rejuvenation etc.

Advantage

- **Diode laser handle**

Painless hair removal, Permanent hair removal, Suitable for all skin tones.

- **Elight handle**

Elight for hair removal ,skin rejuvenation,vascular therapy ,pigmentation therapy,acne therapy.

- **ND YAG Laser handle**

ND Yag laser for tattoo removal, eyeline remova, eyebrow ,embroidery removal. freckle removal.

- **RF handle**

Skin tightening/lifting, wrinkle removal on bod



DIODE LASER

755/808/1064nm diode laser, for all type skin hair removal, painless

755 NM

The Power of Diversity



755nm

808 NM

Up To Speed



808nm

1064 NM

Dedicated to Darker Skin



1064nm





One minute cooling

Suitable for hair removal on any part of the body



Arms



Legs



Face



Armpits



Back




Chest

How does it work

E-light IPL for Hair removal/Skin rejuvenation/Vascular Job/Pigmentation treatment/Acne clearance

Face and body with big and small treatment head.

Through the scientific design of diode laser handle, the high-energy laser beam instantly penetrates the skin and reaches the hair follicle during treatment, output light is absorbed by melanin in the hair follicle, so destroying the hair follicle and reach permanent hair removal. The laser pulse width is set shorter than the time that human body's feeling time, a single shot of light is over, but the human body fell painless, then could achieve the painless treatment feeling, to improve the customer's comfortable experience



Refuse To Repeat And Leave No Scum



1 week-2 weeks
(2 times a week)

75% hair loss
Start to grow slowly
Smaller pores

3 week-4 weeks
(2 times a week)

Hair becomes thinner and
less
Rejuvenate tender
skin
Smooth body

5 week-8 weeks
(2 times a week)

Select according to hair
growth status. Use it
selectively, say goodbye
Repeated hair removal.

IPL Handle

© **High-quality IPL**

Handpiece with super cooling system

© **Optional Spot Size**

15*50mm,15*35mm,12*30mm

© **Long shots time**

IPL:150,000 times

© **7 filters**

420/510/560/585/640/690/755nm

Standard:420/510/560/690nm



420nm

Acne removal
Vascular
treatment

510nm

Pigmentation
therapy,Remove
freckles

560nm

Skin
rejuvenation
wrinkle removal

585nm

Remove Red
blood silk

640nm

Hair removal
(light skin)

690nm

Hair removal
(brown skin)

755nm

Hair removal
(black skin)





ND YAG

5 Treatment Heads

532: Remove colored tattoos.
 755 honeycomb head: used to remove age spots, freckles, and chloasma. little damage to skin
 Adjustable parameter 1064: for tattoo removal, eyebrow clearing
 Non-adjustable parameter 1064: used for removing tattoos and clearing eyebrows
 1320: large area treatment, large spot, fast treatment speed and good effect




ND YAG ADUANTAGE

Imported crystal bonding laser rod, best quality , shot more than 1.5 million times
 Counter on the handle, more easy to get the total shot time
 LED red indicator, more easy to aim on the treatment area.
 5 treatment tips for multi fuctions: tattoo removal, skin rejuvenation, pigment removal, carbon peeling...



RF Handle

Suitable to 6 kinds of therapy area, to make sure the treatment on the corner of eye, and other area on face.



36mm big tripolar rf head


28mm small tripolar rf head

12mm eye-rf

38mm bipolar-rf

28mm bipolar-rf

14mm bipolar-rf



Multifunctional elight ipl diode laser 808nm hair removal nd yag tatoo removal

14inch colour touch screen, man and woman skin, 6 language to choice.

Your language and logo can be added freely.

14 inch true color touch screen

Professional Screen, Easy operation



Delivery time: 3-5 working day after receiving payment.

Payment terms: T/T (bank transfer), Western Union

Warranty:2 years for the whole machine, 1 year for handle

Shipping Way: Door to door service by air, DHL/UPS/TNT

Package: Aluminium box with hard foam

Packing

safe and thoughtful packaging service



Built-in prevention Collision bubble

Two-year warranty for machine handle as consumable,
one-year warranty.

OEM&ODM service for distributor.

Tight packaging

Training: (Video+Manual+Online Service)
Can Keep You Operate It Easily.



3-5 business days delivery

Professional guide:

24 Hours Professional online customer service,