



DIODE LASER HAIR REMOVAL

PAINLESS HAIR REMOVAL, NO SIDE EFFECT

MORE

How Does It Work?

Through the scientific design of diode laser handle, the high-energy laser beam instantly penetrates the skin and reaches the hair follicle during treatment, output light is absorbed by melanin in the hair follicle, so destroying the hair follicle and reach permanent hair removal.

The laser pulse width is set shorter than the time that human body's feeling time , a single shot of light is over, but the human body fell painless, then could achieve the painless treatment feeling, to improve the customer's comfortable experience. Most effective laser hair removal machine.



Three wavelength

755nm 808nm 1064nm diode hair removal technologies in one handle

Three wavelength

755nm 808nm 1064nm diode hair removal technologies in one handle



Depth Of Penetration

Professional Aesthetic Machines Ice Cooling Diode laser Hair Removal Laser Diode Láser Diodo 755 808 1064 Hair Removal

Spot size 12x36mm / 12x24mm Spot Size, for improved absorption of energy at the target follicle.

The only laser with unique and dedicated technology for pain reduction, enhanced efficacy, safety.



■ spot size:12*36



■ spot size:12*24



spot size:12*36



spot size:12*24



spot size:12*36



spot size:12*24

Versatile



The MLKJ laser includes two optimal wavelength 800 nm ice laser diode hair removal systems on one versatile platform.

The MLKJ laser heads enables permanent hair reduction from all areas of the body. Its advanced technology reduces pain while increasing efficacy, enabling fast and comfortable treatments for patients and practitioners.

Efficacy

Perment hair removal for all skin type

Very fair	755nm suitable for light skin hair removal
Fair	
Medium	808nm suitable for neutral skin hair removal
Oliver	
Brown	1064nm suitable for black skin hair removal
Black	




Safety

Almost no skin scattering, no harm to the skin and sweat glands, no scar, no any side-effect


Refuse to repeat and leave No scum



1 week - 2 weeks
2 times a week
- 75 hair loss - Start to grow slowly - Smaller pores



3 week - 4 weeks
2 times a week
- Hair becomes thinner and less rejuvenate tender skin smooth body



5 week - 8 weeks
2 times a week
- Select according to hair growth status use it selectively, say goodbye repeated hair removal

Decreased pain

Patient comfortable by continuously cooling the skin, pain and discomfort are significantly reduced.

Diode laser Hair removal before and after

Hair removal principle

Make Your Beauty Business More Easier



The laser selectively acts on the melanin in the hair follicle, destroying the hair-growing area

Natural hair loss, prevent regeneration, and achieve the purpose of hair removal

Hair removal also stimulates collagen regeneration, shrinks pores and makes skin firm and smooth



Treatment Effect



14-inch 4K screen, smart touch experience

- * Customizable treatment presets
- * With male and female , 6 skin styles to choose from
- * Easy to operateproviding operation skills and teaching services
- * Free logo and language can be customized on the screen

New upgrade system

The newscreen is clear easy to operate.



Our Company & Panking

MLKJ

CUSTOMER SERVICE

15 YEARS FACTORY, 2 YEARS WARRANTY, 102 COUNTRIES, 60,000+ CUSTOMERS.

MLKJ—MAKE YOUR BEAUTY BUSINESS EASIER



Packing

safe and thoughtful packaging service



Built-in prevention Collision bubble

Two-year warranty for machine handle as consumable, one-year warranty.

OEM&ODM service for distributor.

Tight packaging

Training: (Video+Manual+Online Service)
Can Keep You Operate It Easily.



3-5 business days delivery

Professional guide:
24 Hours Professional online customer service,

Contact us to get more information for beauty machines.