

Picosecond

New upgrade

Powerful

Fastly

Safety



Picosecond laser tattoo removal 532 755 1320 1064nm

q switch nd yag laser tattoo removal machine

Product display



Craftsman quality, Reflect the extreme details Rational design

1. The product has a novel appearance and is mainly sold as a hot model
2. Adjustable handle bracket, strong and durable
3. Non-slip handle, ergonomic design
4. Extended treatment connecting tube

**Smallest size**

1. Adopt core components of medical equipment. The size of the machine is smaller.
2. Lower cross-border freight than ordinary machines.
3. The handle is smaller and more convenient to use

5 TREATMENT HEADS

A solution to a variety of skin problems

532: Remove colored tattoos.

755 honeycomb head: used to remove age spots, freckles, and chloasma. little damage to skin

Adjustable parameter 1064: for tattoo removal, eyebrow clearing

Non-adjustable parameter 1064: used for removing tattoos and clearing eyebrows

1320: large area treatment, large spot, fast treatment speed and good effect




5 treatment heads, detachable and rotatable outer wire design

1. 532: Remove colored tattoos.
2. 755 honeycomb head: used to remove age spots, freckles, and chloasma. little damage to skin
3. Adjustable parameter 1064: for tattoo removal, eyebrow clearing
4. Non-adjustable parameter 1064: used for removing tattoos and clearing eyebrows
5. 1320: large area treatment, large spot, fast treatment speed and good effect



Working Principle



TREATMENT HEADS



532: Remove colored tattoos.



755 honeycomb head: used to remove age spots, freckles, and chloasma. little damage to skin



Adjustable parameter 1064: for tattoo removal, eyebrow clearing



1064: used for removing tattoos and clearing eyebrows



1320: large area treatment, large spot, fast treatment speed and good effect

The maximum energy, the best effect


1. Laboratory data shows that the laser rod reaches 2000mj.
2. Precise calculation of laser value
3. The aiming distance of the indicator light is 6-10cm

1064nm

Precise Calculation Laser Value

The aiming distance of the indicator light is 6-10cm

6-10cm



You think you can. you can
No one can stop me

The best laser wand in the industry



Lifespan
1 million times



Imported stone Q-switch cassette mechanism, the entire solid laser, without replacing Q-switch.

Tattoo removal. Remove black, blue, brown and red tattoo pigments on eyebrow, eye line, lip line and other parts of body Remove birth mark, otas nevus and dermal speckle Instantaneous blasting high laser technology, international standard of production, technical test strictly.



Treatment Effect

a solution to a variety of skin problems

	
BEFORE	AFTER
	
BEFORE	AFTER
	
BEFORE	AFTER





Original pigment Traditional laser Picosecond Laser

Before & After

					
Before	Before	Before	Before	Before	Before
					
After	After	After	After	After	After

Screen Display

Super cool touch Experience

6 language for your choice,Special language can be added

OEM service are accepted,such as adjust parameter,include energy/fre/time scope

12.1 inch true color big touch screen

12.1 inch true color touch screen

The display screen is clear at a glance and easy to operate



Technology	Pico second laser / Nd yag laser
Function	Tattoo freckle acne removal / Blackface doll treatment
Handles	5 Heads
Wavelength	1064nm 532nm 1320nm 755nm
Treatment area	Face Body Face Neck Arms Legs....
Frequency	1-10HZ
Voltage	110V/220V 50-60Hz

16 years manufacturer with rich experience in producing beauty machines! 15 years Alibaba golden supplier CE certificates approved TUV certificates approved Personal custom made & OEM

Experienced manufacturer You will cooperate with a factory with more than 15 years of experience in manufacturing beauty equipment, integrating R&D, design, production and sales. We exported to 126+ countries, supporting 60,000+ users, more reliable and more professional.

Professional R&D and newest technology Relying on our professional R&D and technical team, products will continue to be upgraded, you will get the newest machines with latest technologies, and always be at the forefront of the industry, making you more competitive in the beauty industry.

Contact us now to learn more!