


Picosecond laser tattoo removal 532 755 1320 1064nm q switch nd yag laser tattoo removal machine



Craftsman quality, Reflect the extreme details Rational design 1. The product has a novel appearance and is mainly sold as a hot model 2. Adjustable handle bracket, strong and durable 3. Non-slip handle, ergonomic design 4. Extended treatment connecting tube



Smallest size 1. Adopt core components of medical equipment. The size of the machine is smaller. 2. Lower cross-border freight than ordinary machines. 3. The handle is smaller and more convenient to use




5 treatment heads, detachable and rotatable outer wire design


- 1. 532: Remove colored tattoos.
- 2. 755 honeycomb head: used to remove age spots, freckles, and chloasma. little damage to skin
- 3. Adjustable parameter 1064: for tattoo removal, eyebrow clearing
- 4. Non-adjustable parameter 1064: used for removing tattoos and clearing eyebrows
- 5. 1320: large area treatment, large spot, fast treatment speed and good effect



Imported stone Q-switch cassette mechanism, the entire solid laser, without replacing Q-switch.

Tattoo removal. Remove black, blue, brown and red tattoo pigments on eyebrow, eye line, lip line and other parts of body Remove birth mark, otas nevus and dermal speckle Instantaneous blasting high laser technology, international standard of production, technical test strictly. A solutionto avariety of skin problems.





Professional nd yag 2 in 1 laser hair removal machine Multifunctional beauty equipment

- \* USA or Germany Dilas imported laser bar longer life.
- \* More than 50 million shots , more than 3 years don't need to change new handle.
- \* 3 waves 1064nm +808nm + 755nm, 3 in 1 all skin type can be used



Safety Use On All Skin Type 1.Alex 755nm for light skin 2.Speed Diode 808nm for yellow skin 3.Long Pulse 1064nm for darker skin



Screen Display
12.1 inch true color big touch screen
Super cool touch Experience
6 language for your choice, Special language can be added
OEM service are accepted, such as adjust parameter, include energy/fre/time scope



More product details and application methods, welcome to consult us! 7\*24 hours online service!