

[CE Two years warranty diode laser 755 808 1064 diode laser hair removal machine](#)

The screenshot shows the FDA's 510(k) Premarket Notification page. At the top, there is the FDA logo and navigation links for Home, Food, Drugs, Medical Devices, Radiation-Emitting Products, Vaccines, Blood & Biologics, Animal & Veterinary, Cosmetics, and Tobacco Products. A search bar is located in the top right corner. The main heading is "510(k) Premarket Notification" with a sub-heading "FDA Home Medical Devices Databases". Below this, there is a "CDRH SuperSearch" logo and a list of navigation links: 510(k), DeNovo, Registration & Listing, Adverse Events, Recalls, PMA, HDE, Classification, Standards, CFR Title 21, Radiation-Emitting Products, X-Ray Assembler, Medsun Reports, CLIA, and TPLC. The main content area is titled "New Search" and "Back To Search Results". It contains a table of details for a specific 510(k) notification.

Device Classification Name	Powered Laser Surgical Instrument
510(K) Number	K200525
Device Name	Medical Diode Laser Hair Removal System
Applicant	Weifang Mingliang Electronics CO., LTD No.13426, Yuqing Street, High-Tech District Weifang, CN 261000
Applicant Contact	Li Mingliang
Correspondent	Beijing Believe-Med Technology Service Co., Ltd. Rm. 912, Building#15, XiYueHui, No.5, YiHe North Rd., FangShan District Beijing, CN 102401
Correspondent Contact	Ray Wang
Regulation Number	878.4810
Classification Product Code	GEX
Date Received	03/02/2020
Decision Date	05/27/2020
Decision	Substantially Equivalent (SESE)
Regulation Medical Specialty	General & Plastic Surgery
510k Review Panel	General & Plastic Surgery
Type	Traditional
Reviewed By Third Party	No
Combination Product	No

Page Last Updated: 06/01/2020

Has passed Medical CE, ISO13485, FDA seriously

2 New diode laser 755nm 808nm 1064nm depilazion 3 wavelength diode hair removal alma soprano ice platinum

Product Description



Stronger, Faster Two Handles

With 14-inch HD 4K Screen

Option to install two handles for alma soprano ice platinum

Handle cooling



Large Spot Size

Quick and Comfortable

With large spot size of 1.2*3.5cm², for faster treatments and better patient experience.



TEC Plus

Continuous Contact Cooling

TEC is an advanced technology that continuously cools skin and facilitates range-controlled temperature throughout the entire treatment. Patients can now enjoy a completely pain-free, cool and comfortable experience throughout the session.

Features of Alma Soprano Platinum Diode laser hair removal:


- »»Safety:No skin scald, no harm to the skin and sweat glands, no scar, no any side-effect.
- »»Quick: big spot size 12*35 mm², effective and quick treatment on big treatment area
- »»Permanent hair removal: remove all hair on all skin types, including tanned skin.
- »»Crystal temperature: ~5-3°C, reduce the painful and increase comfortable sense in treatment.

Handle features of Alma Soprano Platinum Diode laser hair removal :



- »»**USA laser bar**, 12 bars, high power 1200w
- »»**TEC condenser** cooling system: 2 big 1200W TEC condensers, keep the water constantly at 24°C, ensures the machine working continuously without stopping
- »»**Super Cooling effect**: Sapphire Temperature-5~3°C. Effectively reduce the painful in the treatment and make
- »»**Italy Water Pump**: more better for laser water circle, take away heat, longer life time.
- »»**Two water filters**:two food grade exchange resin,can ensure the water quality for circulation

Different spot size Handle color optional








▶ 10bar




▶ 12bar



▶ 16bar






755 nm for white skin (fine, golden hair removal);



808 nm for yellow neutral skin (brown hair removal);



1064 nm for black skin (black hair removal)

Good Effect of Alma Soprano Platinum Diode laser hair removal

Treatment effect



Screen features of Alma Soprano Platinum Diode laser hair removal

- »»**Touch Screen:** 12 inch easy to operate. 6 language for your choice. Logo customization.
- »»**Color for choose,** also your LOGO could be put on screen and machine shell for free
- »»**Metal inner structure,** more stable, and suitable for long distance delivery.

Screen image

12 inch touch screen easy to operate.



Model	ML-808 V19 Alma Soprano Platinum Diode laser hair removal
Laser type	808nm /755 808 1064
laser bar	12bar , 1200w, imported from USA
Display system	12 inch true color touch screen
Laser handle shots	40 million times
Energy	10~150 j/cm2, adjustable
Pulse Duration	40~800ms adjustable
Frequency	1~10Hz
Cooling system	TEC condenser
	Air and close water circulation
	Treatment head with semiconductor
	cooling temprature adjustable
Power	3000w
Spot size	12*35mm
Crystal temperature	-10~3°C
Power requirement	AC220/110V, 50/60HZ
N.W.	65KG
G.W.	80KG
Box size	122*65*68cm

Alma soprano ice / Alma Titanium laser / Alma soprano ice platinum diode laser for hair removal

Company advantage

Service users all over the world



Packing

safe and thoughtful packaging service



Built-in prevention Collision bubble

Two-year warranty for machine handle as consumable,
one-year warranty.

OEM&ODM service for distributor.

Tight packaging

Training: (Video+Manual+Online Service)
Can Keep You Operate It Easily.



3-5 business days delivery

Professional guide:
24 Hours Professional online customer service,



FAQ

Q1. Is it a reliable company?

MLKJ reply: Weifang Mingliang Electronics Co.,Ltd. Is established in 2005,engaged in R&D,sales and after-sales services of aesthetic devices & medical laser equipments,even part of IT area and real estate sector.We have a professional team integrating with optics, machinery,electricity and medicine,to keep us ahead on in the field. Our factory covers 5000m² and has more than 100 employees.

Q2.What about the delivery?

MLKJ reply:we can offer the door to door services by air or by sea, just depend on your actual request.Such as DHL,UPS,TNT, FEDEX...By air ,and the sea transportation.

Q3.What's the delievery time?

MLKJ reply:3 working days.Huge stock.

Q4.What's package?

MLKJ reply:strong and beautiful aluminum alloy case/carton case/wooden case

Q5:Do you have any timely technology supports?

MLKJ reply:We have a professional technology supporting team for your timely services.We prepare the technical documents for you, also you can contact us by telephone,webcam,online chat(googletalk,skype,yahoo,whatsapp,viber...)

Q6:What's the payment method?

MLKJ reply:T/T,WESTERN UNIN,ESCROW

Q7.If the machines broken during the shipment,will you support us?

MLKJ reply:Generally,our customer will purchase the insurance before shipment,avoid any money lost from rough transportation. We will help you to solve the claims with the insurance companies to collect your money back in time.



Contact us

If you want more detail,please contact me Jennifer Dong.

You can Click **Send inquiry** ,convenient for you to get more information quickly.

2 Years Warranty

- ✓ Free training
- ✓ User manual
- ✓ Installation Precautions
- ✓ Treatment parameter setting

