An advertisement for the GENES LASER hair removal system. On the left, a woman with curly hair is smiling. On the right, a black laser machine on wheels is shown with a '50% OFF' sticker. The background is dark with orange geometric shapes.

GENES LASER

**NEWEST COMING, FASTEST
STRONGEST, DOUBLE HAIR
REMOVAL SYSTEM**

1200W + 1800W TITANIUM COMING

[VIEW MORE](#)

*1800w High power diode laser with nice cooling technology, for hair removal with painless!

*Diode laser enable the light to penetrate deeper into skin and safer than other laser. because it can avoid the melanin pigment in the skin's epidermis, we can use it for the permanent hair removal of all color hairs on all 6 kinds skin types.

*Suitable for any unwanted hair on areas like face, arms, armpits, chest, back, bikini, legs...It also have the skin renewing.

*Frequency 1-10hz. treatment fastly!!! Machine for fast and permanent hair removal . Painless!!

Handle display

- ① Optional 12*24mm/12*36mm light spot size.

Different spot size meet different body parts hair removal needs .



*spot size : 12*24mm*



*spot size : 12*36mm*

- ② 10bars 1200w / 15bars 1800w. 10million—40million shots.

Effective for permanent hair removal

- ③ 755nm and 808nm and 1064nm , suits for all skin and hair types.

Cooling effect

The handle has a cooling effect that improve treatment comfort.



*spot size : 12*24mm*



*spot size : 12*36mm*

Waves laser combines 755nm 808nm 1064nm, suits for all skin and hair types treatment




755nm for white skin (fine/ golden hair) hair removal


808nm for yellow (natural) skin hair removal

1064nm for black skin (black hair) hair removal


diode laser hair removal machine

working principle


-  755 nm for white skin (fine, golden hair removal);
-  808nm for yellow neutral skin (brown hair removal);
-  1064nm for black skin (black hair removal)



A natural




B energy absorbed by hair



C hair follicle died

working principle


Triple wavelength laser spots work together




Before and after effect comparison

Scope of application

Suitable for any unwanted hair on areas like face, arms, armpits, chest, back, bikini, legs...




CE and TUV CERTIFICATE


**U.S. FOOD & DRUG
ADMINISTRATION**

[Home](#)
[Food](#)
[Drugs](#)
[Medical Devices](#)
[Radiation-Emitting Products](#)
[Vaccines & Blood](#)

510(k) Premarket Notification

[FDA Home](#)
[Medical Devices](#)
[Databases](#)



510(k) | Defenses | Registration & Listing | Adverse Events | Recall
CFR Title 21 | Radiation-Emitting Products | X-Ray Access

New Search

Device Classification Name	Crossed Laser Gun
510(k) Number	K090205
Device Name	Medical Grade Laser
Applicant	Wavelength Medical 140 N. 14th, Youngstown, OH 44460
Applicant Contact	L. Mingling
Correspondent	Beijing Bellevue Medical Room 512, Building 1 Fangshan District Beijing, CN 102404
Correspondent Contact	May Wang
Registration Number	524, 4613
Classification Product Code	524
Date Received	03/02/2003
Decision Date	05/27/2003
Decision	Substantially Equivalent
Regulation Medical Specialty	General & Plastic Surgery
510k Review Panel	General & Plastic Surgery
Type	Traditional
Reviewed By Third Party	No
Combination Product	No


Measurement Systems Certification Body No. 2002
Institute for Testing and Certification
T. Bělá 206, 754 21 Zlín, Czech Republic
www.itc.cz

CERTIFICATE
No. 10 0144 SJ

We confirm on the basis of a performed audit that company

Weifang Mingliang Electronics Co.,Ltd.
338Wang Fang Dongjia Street Wafang City
Company Reg. No.: 3707000000000000

has implemented and documented a functional quality management system
in compliance with the requirements of the standard

EN ISO 13485:2003

Covering the following activities:

**Design & Development, Production and Distribution of IP, Phototherapy System,
Model: VAG-G Softlaser Laser and RF Device**

The Certificate is based on the basis of the results mentioned in Audit Report No.
2004000000000000. The Certificate validity is conditioned by positive results of surveillance audits
which the certified company undertakes in regular

Paul Vög

By: Paul Vög
Head of Certification Body


Date of Issue: 01.01.2004
Valid until: 31.01.2004
Date of the last certification audit: 04.07.2007



SHIPPING

Welcome to visit our factory



In terms of following trade details:

- 1)The shipping cost are based on DHL cost to your country
- 2)Delivery time:3-5 days after receiving payments.
- 3)Payment terms:T/T(bank transfer),Western Union.
- 4)Warranty:2 years for the whole machine
- 5)Warranty Policy:one month guarantee replacement (we bear all shipping cost),2years warranty,lifelong maintenance.
- 6)Packing method: Aluminum box with hard plastic foam.
- 7)Free-training to customers:Free We could offer DVD and User's manual and 10 Hours online service to our customers
- 8)Certifications:CE TUV ISO 13485
- 9)MOQ:1 UNIT
- 10)The quotation here above is valid within 1 month