





- \*1800w High power diode laser with nice cooling technology, for hair removal with painless!
- \* Diode laser enable the light to penetrate deeper into skin and safer than other laser. because it can avoid the melanin pigment into the skin's epidermis, we can use it for the permanent hair removal of all color hairs on all 6 kinds skin types.
- \* suitable for any unwanted hair on areas like face, arms, armpits, chest, back, bikini, legs...It also have the skin renewing.
- \* Frequency 1-10hz. treatment fastly!!!

Machine for fast and permanent hair removal . Painless!!



Product Paramenters



Model	808V19
Laser Type	808nm diode laser or 755/808/1064nm 3 waves combined
Display	14 inch 4K touchable screen
Wavelength	755/808/1064
shots	20-40 million shots
Energy	10~150 j/cm2, adjustable
Frequency	1-10HZ
Cooling	Compressor+air+water
Power	3500W
Spot size	12*24mm or 12*36mm for choose
Voltage	110V/ 60HZ or 220V/50HZ for choose
Net Weight	95Kg



Gross Weight 105Kg  
Machine size 53\*40\*105cm  
Box size 125\*56\*70.4  
Handle display

## Handle display

① Optional 12\*24mm/12\*36mm light spot size.

Different spot size meet different body parts hair removal needs .



*spot size : 12\*24mm*



*spot size : 12\*36mm*

② 10bars 1200w / 15bars 1800w. 10million—40million shots.  
Effective for permanent hair removal

③ 755nm and 808nm and 1064nm , suits for all skin and hair types.

## Cooling effect

The handle has a cooling effect that  
improve treatment comfort.



*spot size : 12\*24mm*



*spot size : 12\*36mm*

4cm<sup>2</sup>



## Three Simultaneous Wavelengths in One Applicator

Plus Exceptionally 4cm<sup>2</sup> Large Sport Size



The lowest handle  
temperature can reach -27°C



# For all skin types

**Triple wavelength** laser work together

755  
nm

808  
nm

1064  
nm

755nm/808nm/1064nm

755nm  
808nm  
1064nm

**Waves laser combines 755nm 808nm 1064nm, suits for all skin and hair types treatment.**

755nm for white skin (fine/ golden hair) hair removal.

808nm for yellow (natural) skin hair removal.

1064nm for black skin (black hair) hair removal.

Before and After ice platinum diode laser hair removal machine




# working principle



755 nm for white skin (fine, golden hair removal );



808nm for yellow neutral skin (brown hair removal);




1064nm for black skin (black hair removal)



**A** natural




**B** energy absorbed by hair



**C** hair follicle died

# working principle

Triple wavelength laser spots work together




Before and after effect comparison

## Scope of application

Suitable for any unwanted hair on areas like face, arms, armpits, chest, back, bikini, legs...













**We help you**  
Goodbye to the embarrassment



The sign displays a 3x2 grid of 'Before' and 'After' photos showing hair removal results on arms, legs, and chest.

Before	After
	
	
	

Screen system display

# System Function



- 4k high definition screen system
- Different parts can be selected
- 6 skin type Any language could be added
- Logo could be added for free
- Easy to operate
- Stable operation

## The screen of Alma soprano titanium

- 14 inch big screen, quick touch, easy operate
- Male & female, 6 body parts, 6 skin types
- 6 language options, your native language can be added

## Our certificate

All the process of the producing is under ISO 13485 Quality system and match with CE certification request.



BV-France's most authoritative  
certification body BV Group certification



TUV-Germany's most authoritative  
certification body Rheinland certification



SGS-Switzerland's most authoritative  
certification body SGS certification

