

# Diode Laser

ice Hair removal



1200 1800

1. .  
24 \* 12 : .2mm 12 \* 36mm  
10-1 : .3Hz  
.  
.  
.  
.



□□□□□□□ □□□□□ □□□□ □□□□ □□□□ □□□□

□□□□□□ □□ □□□□□ □□□□□ □□□ □□□□ □□ □ □□□□□ □□ □□□□ □□□□ □□ □□ 24 \* 1.12 □□ □□□□□□ □□□□□□ □□□□□□□□

□□□□□□□□ □□□□□□□□ □□□□ □□□ □ □□□□□ □□ □□□□□ □□□□ □□ □□ 36 \* 2.12 □□□□□□□ □□□□□□□□

□□□□□□□ **1064** □□□□□□□ **808** □□□□□□□ **755** :□□□□□□□ □□□□□□□ □□□□□ □□□□□

755 .□□□□□□□ □□□□□□□ □□□ □□□□□□ □□□□□□□□  
 .□□□□□□□□□ □□□□□ □□□ □□□□□□ □□□□□□□□ 808  
 1064nm □□□□□□□□ □□□□□□ □□□ □□□□□□ □□□□□□□.

Medical 1200W/1800W

50 Million Shots

USA / Germany Imported Laser Bar



755nm+808nm+1064nm



ALEX 755nm

**For white skin**

(fine, golden hair removal )



SPEED 810nm

**For yellow neutral skin**

(brown hair removal)



YAG 1064nm

**For black skin**


(black hair removal)

10 1800 15 1200:

\* ( ) .

50 20 :

1064 755808 \*nm



000 0000

1). 0000 14 00000 000000 0000 0000 :000000 00000

.0000000 000 0 000000 000000 00 000000 000000 0000 :0000000 0000 (2

00000 (3 OEM: 0000000 00000 000000 0 0000000 00000000 00000 0000 000000 000000 00000.