



□□□□□□□□ □□□□□□ □□□□ □□□□ □□□□ □□□□

.□□□□ □□□□ □□□□□□ □□□□ □□□ □□□ □□□ □□ □□□□ □□ □□□□ □□□□ □□□ □□ □□ 24 * 1.12 □□ □□□□□□□□ □□□□□ □□□□□□□□ □□□□□□□□ □□□□□ □□□ □□□ □□□ □□□□□ □□□□□ □□□ □□ □□ 36 * 2.12 □□□□□□□□ □□□□□...

Arm hair removal



Face hair removal



Lip hair removal



Back hair removal



Leg hair removal



Bikini hair removal



Underarm hair removal



1064nm 808nm 755nm : 1064nm 808nm 755nm

755nm . 1064nm 808nm 755nm . 1064nm 808nm 755nm . 1064nm 808nm 755nm .

Medical 1200W/1800W

50 Million Shots

USA / Germany Imported Laser Bar



755nm+808nm+1064nm



ALEX 755nm

For white skin

(fine, golden hair removal)



SPEED 810nm

For yellow neutral skin

(brown hair removal)



YAG 1064nm


For black skin

(black hair removal)

10 1800 15 1200 :

* () . 50 20 :

1064 755808 *nm



000 0000

1). 0000 14 00000 000000 0000 0000 :000000 00000

.0000000 000 0 000000 000000 00 000000 000000 0000 :0000000 0000 (2

00000 (3 OEM: 0000000 00000 000000 0 0000000 00000000 00000 0000 000000 000000 00000.