

Diode Laser

ice Hair removal



1200 000000 0000000 000000 0000 1800 0000

1. 000000 000000 00000 0000 000000 0000 0000 0000
24 * 12 0000 0000 :000000 .2mm 12 * 36mm 000000
10-1 0000 :00000000 00000000 .3Hz 0 00000000 00 00000000 0000 000000.
00000000 0000





上海圣兰西医疗美容有限公司
上海圣兰西医疗美容有限公司

Arm hair removal



Face hair removal



Lip hair removal



Leg hair removal



Back hair removal



Underarm hair removal



Bikini hair removal





□□□□□□ **1064** □□□□□□ **808** □□□□□□ **755** :□□□□□□ □□□□□□ □□□□ □□□□

755 .0000000 0000000 000 0000000 00000 00000000
 .00000000 0000000 000 0000000 00000 0000000 808
 1064nm 00000000 0000000 000 0000000 000000.

Medical 1200W/1800W

50 Million Shots

USA / Germany Imported Laser Bar



755nm+808nm+1064nm



ALEX 755nm

For white skin

(fine, golden hair removal)



SPEED 810nm

For yellow neutral skin

(brown hair removal)



YAG 1064nm


For black skin

(black hair removal)

10 1800 15 1200 :

* () . 50 20 :

1064 755808 *nm



000 0000

1). 0000 14 00000 00000 0000 0000 :000000 00000

.0000000 000 0 000000 00000 00 00000 000000 0000 :000000 0000 (2

0000 (3 OEM: 0000000 0000 00000 0 0000000 0000000 0000 000 000000 000000 00000.

4K 14 inch HD screen system

