

NEW ARRIVAL



C2 Vertical Magnetic Thin

- ◆ BUILD MUSCLE
- ◆ BURN FAT

The newest air-cooling system
4 working handles





- 1. 30 5.5 =
- 2.
- 3. 19
- 4. 16
- 5.
- 6.
- 7.
- 8.



Operation



Hamstring



Thighs



Abdomen




Arm



Hip



BEFORE AND AFTER



Before

5 Weeks After

8 Weeks After



Before

After



Before

After







15 Years Laser Beauty Device Manufacturer

Make your beauty business easier



Smayer is a leading manufacturer of laser beauty devices, providing high-quality products and services to beauty businesses worldwide. With over 15 years of experience, Smayer has established a strong reputation for innovation and reliability. Our products are designed to be safe, effective, and easy to use, making them the ideal choice for beauty professionals. We offer a wide range of laser devices, including hair removal, skin rejuvenation, and tattoo removal systems. Our commitment to excellence and customer satisfaction has made Smayer a trusted partner for beauty businesses globally.

Our Factory Show

**15 Years
Professional Factory**




**Professional
Production R & D
Sales Department**



**Mass production
of various styles**



**Real Feedback
of Alibaba**



Certificate display

Certificate is the best guarantee for you



□□□□□□ □□□□□□

Packing

Safe and thoughtful packaging service



Tight packaging

The machine shall be made of waterproof materials and packaged

**Anti collision foam board
inside aluminum box**



3-5 business days delivery

Professional guide:

24 Hours Professional online customer service

Our company has a professional packing team, 7-3 : packing service
We can provide you with professional packing service, DHL, TNT, Fedex, etc.
We can provide you with professional packing service, etc.
We can provide you with professional packing service, etc.

CONTACT US



In terms of following trade details:

- 1)The shipping cost are based on DHL cost to your country
- 2)Delivery time:3-5 days after receiving payments.
- 3)Payment terms:T/T(bank transfer),Western Union.
- 4)Warranty:2 years for the whole machine
- 5)Warranty Policy:one month guarantee replacement (we bear all shipping cost),2years warranty,lifelong maintenance.
- 6)Packing method: Aluminum box with hard plastic foam.
- 7)Free-training to customers:Free We could offer DVD and User's manual and 10 Hours online service to our customers
- 8)Certifications:CE TUV SFDA ISO 13485
- 9)MOQ:1 UNIT
- 10)The quotation here above is valid within 1 month

一、产品概述

1.1 产品名称: 智能血压计

1.2 产品型号: 7 系列智能血压计

1.3 产品规格: 2.0 英寸显示屏

1.4 产品特点: 本产品采用先进的传感器技术, 能够准确测量血压, 并具有数据存储功能。

1.5 产品用途: 适用于家庭、医院、诊所等场所。

1.6 产品优势: 本产品具有操作简单、测量准确、体积小、重量轻等优点。

1.7 产品包装: 本产品采用坚固的塑料盒包装。

1.8 产品价格: 4.0 元/台

1.9 产品售后服务: 本产品提供 2 年质保, 终身维护。

1.10 产品认证: 5.0 认证

1.11 产品运输: 本产品采用 DHL 快递运输, 运费自理。