



755 nm / 808nm / 1064nm


- * .
- . *
- . *
- * *
- * 1064 808 755 808 *

Suitable for hair removal on any part of the body,
Such as arms, legs, face, armpits, back, chest, etc.



• Cooling effect •

The handle is frosted for one minute,
making your treatment more comfortable and painless






Soon

Permanent hair removal

Three Waves 755nm 808nm 1064nm

755nm (Alexandrite) 808nm (Speed) 1064nm (YAG) .






ALEX 755nm
For white skin
(fine, golden hair removal)




SPEED 810nm
For yellow neutral skin
(brown hair removal)



YAG 1064nm
For black skin
(black hair removal)



Treatment Comparison





Emergency stop button

12 inch big screen quick touch easy operate

Diode laser handle

Vents

Plug

1. 2. 3. 4. 5. 6.

1. .
2. .
3. .
4. .
5. .
6. .

!