



755 nm / 808nm / 1064nm

- *.
- *
- *
- *
- 1064 808 755 808 *

Suitable for hair removal on any part of the body,
Such as arms, legs, face, armpits, back, chest, etc.



• Cooling effect •

The handle is frosted for one minute,
making your treatment more comfortable and painless






Soon

Permanent hair removal

Three Waves 755nm 808nm 1064nm

755nm (Alexandrite) / 808nm (Speed) / 1064nm (YAG) .
 808nm (Speed) .
 1064nm (YAG) .






ALEX 755nm
For white skin
(fine, golden hair removal)



SPEED 810nm
For yellow neutral skin
(brown hair removal)



YAG 1064nm
For black skin
(black hair removal)



Treatment Comparison

