

DIODE LASER



ICELLO



MLKJ . **MLKJ** .

MLKJ . **MLKJ** .

□ FDA & **MLKJ** CE

□ 3 **MLKJ** - 1
 □ Nd: YAG lasers - **MLKJ**
 □ **MLKJ**

□□□□ □□□□□□□□ □□□ □□□□□□□□ □
□□□□ □□□ □□□□ □□□□□ 50□!

□□□□□□ □□□



- 1-□□□□□□□ □□□□□□ □□□□□□ □□□□ □□□ □ □□□□□ □□□□ 32- □□□ □□□□□□ □□□ □□□□□ :□□□□□□ □□ □□□□□□
- 2- □□□□□ □□□□□ □□□□ □□ □□□□ □□□□□□ □□□□□□ 3-2 :□□□□□□□ □□□□□□□ □□□□□□
3. □□□□□□ □ □□□□□ □□□□□□ □□ □□□□□□□ □□ □□□□□□ □□ □□□□□ □□ □□□□□ □□□□□□
- 4- □□□ □□ □□□ □ □□□□ □□□□□ □□□□□□ □□□□ □□□□ □□□□□□□ □□□□□□ □□□□□□ :□□□□□ □□□ □□□ □□□□□□ □□□□□□ □□□□□□□ □□□□□□□□□□□□.



755/808 / 1064nm .□□□□□□□□ □□□□□ □□□□□ □□□ □□□□
808<< single wave or 755nm/808nm/1064nm perfect triple wavelengths
>>America Imported Emitter, strong energy
>>10/12 pieces or 16 pieces Germany Dilas original laser bars,1200W/1800W power is capable to meet all demands of hair removal on the whole body


Medical diode laser

50000000+shots
can treat 10000+clients

- *Not easy to be burned
- *Low maintenance cost
- *Stable quality



- >>1-10Hz fast frequency,faster, smooth and effective treatment
- >>Super sapphire touch cooling system ensures painless and comfortable treatment
- >> Italy imported stainless water pump, it can keep machine working long time without stop. Make sure the heat were removed very quickly
- >>Double filters to keep water quality



The real experience is visible

Through live test verification ,safe ,long,painless
and effective



4 0000 0000 14 0000K 0 0000 000 00000

* 6 00000 00 000000000 0000 000000

00000000 00000 / 00000000 0000000 0000 000 00000000 000000 *

00000000 000 00000000 00000000 00000000 000000 00000 *



You'll get

Two major service guarantees
Worry-free after-sales

01 Professional after-sales group

24 hours offers the solution
Solve problems quickly with your active cooperation

02 Free maintenance

Warranty period, Send your replacement
parts for free

本公司产品符合 CE 认证，符合 ISO13485 标准，符合 FDA 认证。15 年质保，终身维护，让您无后顾之忧！

Establishment Registration & Device Listing

FDA Home Medical Devices Databases



New Search Back To Search Results

Establishment:
 WEIFANG MINGLIANG ELECTRONICS CO.,LTD.
 Floor 3, Building 4,Second Accelerator, Yuqing Eastern Street, High-Tech District
 Weifang Shandong, CN 261000
Registration Number: 3007549439
FEI Number*: 3007549439
Status: Active
Date Of Registration Status: 2022

Owner/Operator:
[Weifang Mingliang Electronics CO.,LTD.](#)
 Floor 3, Building 4,Second Accelerator, Yuqing Eastern Street, High-Tech District
 Weifang, Shandong CN 261000
Owner/Operator Number: 10082104

Official Correspondent:
 Pingping Lv
 Floor 3, Building 4,Second Accelerator, Yuqing Eastern Street, High-Tech District
 Weifang, Shandong CN 261000
Phone: 86-10-56335780

US Agent:
 Wen Xu
 4611 Maydew Mews Way Ellicott City MD
 Ellicott, MD US 21042
Phone: 443 8968583 Ext
Email: Infor@Believe-Med.Com

* Firm Establishment Identifier (FEI) should be used for identification of entities within the imports message set

Page Last Updated: 05/09/2022

Note: If you need help accessing information in different file formats, see Instructions for Downloading Viewers and Players.

Language Assistance Available: Español | 繁體中文 | Tiếng Việt | 한국어 | Tagalog | Русский | العربية | Kreyòl Ayisyen | Français | Polski | Português | Italiano | Deutsch | 日本語 | فارسی | English



Authentication

FDA, NQA13485, CB certification and more Complete certification materials



60000+ Customers

We have professional after sale team and we provide 7*24 services



2 Years Warranty

2 years warranty for each machine, All questions will get reponse within 24 hours

24 x 7 服务支持！