



# Soprano ICE PLATINUM

Medical beauty machine

Medical diode laser hair removal  
For all skin types,  
triple wavelength laser work together

- . 1064nm 808nm 755nm 3-wave laser
- . 1064nm 808nm 755nm 3-wave laser \* (skin types 1-6)
- . 1064nm 808nm 755nm 3-wave laser \*
- . 1064nm 808nm 755nm 3-wave laser \*
- 1064nm 808nm 755nm 3-wave laser \*
- 1064nm 808nm 755nm 3-wave laser 808nm 1064nm 755nm \*

**350g**  
Ultra light handle



12x36mm

**Handle**

1800W diode laser hair removal



12x24 mm


**Handle**

1200W diode laser hair removal

Three Waves 755nm 808nm 1064nm 1064nm 808nm 755nm 3-wave laser 808nm 1064nm 755nm \*

755nm 1064nm 808nm (1064nm / 808nm 755nm) 1064nm 808nm 755nm \*

.( ) 808 1064nm . ( )




**ALEX 755nm**  
**For white skin**  
(fine, golden hair removal )



**SPEED 810nm**  
**For yellow neutral skin**  
(brown hair removal)



**YAG 1064nm**  
**For black skin**  
(black hair removal)



# Treatment Comparison

