

# Hair Removal + ND YAG

2in1 Multifunctional combination machine  
one machine and two effects



2in1 Multifunctional combination machine

2 in 1 Multifunctional combination machine one machine and two effects 1 in 1 nd yag

2in1 Multifunctional combination machine

## 2in1 Multifunctional combination machine



- ✓ Hair removal
- ✓ Toner whitening
- ✓ Mole removal




- ✓ Tattoo removal
- ✓ Freckle
- ✓ Birthmark removal



# Real shot display


Higher-definition screen, stronger cooling effect,  
more effective treatment




1. 4. 1 RF Elight IPL 50 500000 500000 500000 500000  
 . 50 500000 500000 500000 500000 500000 500000 500000  
 \* Elight 30000 150000 500000 500000 500000  
 \* Nd yag: 1064 5 1064 1320 nm  
 \* RF: 6

## *Three wavelengths work together*

**12\*24mm<sup>2</sup>** {



**12\*35mm<sup>2</sup>** {



**Suitable for hair removal on any part of the body,  
Such as arms, legs, face, armpits, back, chest, etc.**







# HANDLE-ADVANTAGE



## specification

- Multiple choice  
fixed 1064nm & 532nm, adjustable 1064nm & 532nm (adjust 1-7 mm) for tattoo removal, nevus removal etc.
- Multifunction  
1320nm for skin rejuvenation
- Indicator light  
Dot position is the indicator light position dot more accurate



Adjustable  
532nm



Adjustable  
1064nm



Black doll  
1320nm



Fixed  
532nm



Fixed  
1064nm



## Before and after



## 14inch true touch screen

6 kinds skin type for choice. accept ODM OEM Service.

The new **14-inch** screen is clear easy to operate.



□□□□□□□□ □□□□□□ □□□□□□ □□□□□□ □□ □□□□□□□□ □□ □□□□□□ □

□□□□ □□ □□□□ □□□□□ □□□□

**24** □□□□□□□ □□□ □□□□ !!