

# Diode Laser

ice Hair removal



**1200** **1800**

1. .  
24 \* 12 : .2mm 12 \* 36mm  
10-1 : .3Hz



1064nm 808nm 755nm  
 1064nm 808nm 755nm

1064nm 808nm 755nm

1064nm 808nm 755nm 24 \* 1.12

1064nm 808nm 755nm 36 \* 2.12

1064nm 808nm 755nm

755nm 808nm 1064nm

Medical 1200W/1800W

**50 Million Shots**

**USA / Germany Imported Laser Bar**



**755nm+808nm+1064nm**



**ALEX 755nm**

**For white skin**

(fine, golden hair removal )



**SPEED 810nm**

**For yellow neutral skin**

(brown hair removal)



**YAG 1064nm**


**For black skin**

(black hair removal)

**10 1800 15 1200 :**

\* ( ) . 50 20 :

**1064 755808 \*nm**



000 0000

1). 0000 14 00000 000000 0000 0000 :000000 00000

.0000000 000 0 000000 000000 00 000000 000000 0000 :0000000 000 (2

0000 (3 OEM: 0000000 00000 000000 0 0000000 00000000 00000 000 000000 000000 00000.