



Medical Diode Laser Hair Removal


1800W+1200W handles



000000 0000



000000 000






Laser show Handle laser bar optional



④ 10bar

④ 12bar

④ 16bar



□□□□□□ □□□□

Suits for all skin and hair types



808nm

1064nm



- **755nm for white skin** (fine, golden hair removal)

● **808nm for yellow neutral skin** (brown hair removal)


- **1064nm for black skin** (black hair removal)



- 1. 1000000 1000000 1000 1000 10 100000 10 10000 15 10 1000 10000 .1000000.**

[illegible][illegible]

Treatment Comparison



□□□□ OEM

1. □□□□□ □□□ □□□□
2. □□□□□□□ □□□□□□□ □□□ □□□□□ □□□□

customer visit



Packing

Safe and thoughtful packaging service



Tight packaging

The machine shall be made of waterproof materials and packaged

**Anti collision foam board
inside aluminum box**



3-5 business days delivery

Professional guide:

24 Hours Professional online customer service