



808+ND

Diode laser hair removal
 ND YAG laser tattoo removal

Nd yag

1. Nd yag laser 1 or 2 sessions
2. Shots per session: 40-100
3. Nd yag: 1064nm / 532nm







□□□□□ □□□

Three wavelengths work together

12*20mm² {



12*35mm² {



Suitable for hair removal on any part of the body,
Such as arms, legs, face, armpits, back, chest, etc.



755 808 1064 波長のレーザー光線が、毛根を破壊し、毛の成長を抑制します。
効果は持続的であり、40-20 回の施術で、毛の再発を抑制します。




Before

After



Before



After



Before



After

5 therapeutic **Heads**

- 532 & 1064: removing tattoo, eyeliner and lip line
- 1320: black-faced doll whitening and whitening skin
- Adjustable 532 & 1064: removing large-area tattoos



Adjustable
532nm



Adjustable
1064nm



Black doll
1320nm



Fixed
532nm




Fixed
1064nm

□□□□ □□□□ □□□□ □□□□:



1064 □ □□□□□ □□□□ 5 □□□□nm / 532nm 1064 □ □□□□□□ □□□□nm / 532nm 1320 □ □□□□nm □□□□.

Working Principle



1. 000000 000000 6 0 000000 000000 0 0000 / 0000.
2. 0000 0000 000000000000 000000.
3. 00000000 00000000 0000 OEM (0000 / 0000...)

System Design



808 Mode:Operation interface

Laser Mode:Tattoo removal

808 nd yag laser

808 nd yag laser Nd Yag 808

2000

SHR 1-10

808 12 * 20 35 * 12

Nd yag 5

10-50 2 /

1-10

-21

12 TFT LCD

220 60 / 110 50 /

116 * 60 * 61

90

90



15 Years Laser Beauty Device Manufacturer


Make your beauty business easier



1. Weifang Sunshinekj Electronics Co., Ltd (Chinese: MLKJ) is a professional manufacturer of laser beauty devices. We have 15 years of R & D experience in the field of laser beauty devices. We are committed to providing high-quality products and services to our customers.

Premium OEM Factories

High-performance manufacturing capacity
R&D capability for customization
Professional certifications and approvals



Packing

Safe and thoughtful packaging service



Tight packaging

The machine shall be made of waterproof materials and packaged

**Anti collision foam board
inside aluminum box**



3-5 business days delivery

Professional guide:

24 Hours Professional online customer service