



40K CAVITATION 5IN1

RF+VACUUM+CAVITATION
+LASER PADS SLIMMING MACHING

START



Advantage^e



Advantage

- **portable**

5 in 1 salon beauty system
saving shipping cost

5 handles

for slimming, skin tightening, wrinkle removal

- **Advantage**

high power cavitation 300W more effective

vacuum head help a lot for supus fat cell

2 years warranty

multi-color for option





Cavitation

Frequency:4Mhz

Power:150W

Burning explode fat cell



Monopolar RF

Frequency:40khz

Power:300W

Body shape fat reduction



Big tripolar RF

Frequency:4Mhz

Power:150W

Skin tightening lifting wrinkle removal



Small Tripolar RF

Frequency:1Mhz

Power:150W

Skin tightening lifting wrinkle removal



Vacuum

Frequency:1MHz

Pressure:-90kpa

Slimming Remove surplus fat cell

□□□

□□□□ TFT LCD

□□□□□□□□

□□□□□□□

□□□□□□□ □□□□

□□□□□□□ □□□

□□□□□□

□□□□□□□ □□□ □□□□ □□□□□□ RF

□□□□□ 8 □□□□ □□□□□□ □□□□ □□□□ □□□□

□□□□ □□□□□□□ RF (RF □□□□□ □□□□□ RF □□□□□□□ □□□□□)

□□□□□□□ □□□ □□□□□□□ □□□□□□□ □□□ □□□ RF

□□□□□□□ □□□ □□□

□□□ RF □□□□□ □□□□□□

□□□ RF □□□□□) □□□□□□□□ □□□□□□□

□□□ RF □□□□□) □□□□□□□□ □□□□□□□

40 □□□□ □□□□□

300 □□□

RF 100000	1-10 000 0000 0000) 000 150 (00000000 0000) 0000 00000)
00000000	
000000000 RF	1-10 000 0000 0000) 000 150 (00000000 0000) 0000 00000)
0000000 000	100 0000000 00000
00000	00 38 * 31 * 47 :0000000000 000
	00 46 * 43 * 54 :000000000000 0000 000000 000
000000000 000000	220 000000 00000 00000 110 / 000000 00000 00000
00.	15 000
00000000.	25 000

Treatment before & after compare




00000000 00000000 000000:

- >Suitable for 6 different color skin,including the dark skin.
- >The sapphire temperature cools down 0~3°C,clients feel comfortable and painless in whole treatment .
- >Long continuous working time,ensure the works in Clinic/salon use.

AMAZING EFFECT

Before & After



Before


Five weeks later

Eight weeks later



Before

After



Before

After



Before

After



Before

After

Free logo service



□□□□ □□ □□□□!