



From Israel Technology

FOR ALL SKIN TYPE HAIR REMOVAL

- 100% Ice Cooling
- 100% Permanent
- 100% Painless

.0000 000000 00000 00000 000
 .000000 00000 00 000 000 00000 .0000 00 *
 .000000 00000 00 000 000 00000 *
 .00000 000 *
 0 00000 000000 000000 000000 0000000 *
 0 0000000 1064 808 755 00 0000000 808 000000 000000 000000 000000 *

350g

Ultra light handle



12×36mm

Handle

1800W diode laser hair removal



12×24 mm


Handle

1200W diode laser hair removal

0000 0000 **Three Waves 755nm 808nm 1064nm** 0 000000 000000 000000 00000 0000 00000 0000.

755 .000000 0000000 (0000000 / 0000000 000000) 00000000 00000000 00000000

.() 808 1064nm () .




ALEX 755nm
For white skin
(fine, golden hair removal)



SPEED 810nm
For yellow neutral skin
(brown hair removal)



YAG 1064nm
For black skin
(black hair removal)



Treatment Comparison

