

Stronger, Faster

×

14-inch HD 4K Screen



000000 000

----- 0 000000 000000 000000 0000 0 000 0000 00 0 000000 000000 0000 000 0000 0000 0000 000000 000000

----- 000000000 000000 000000 000000 0000000000 00000 000 000000 00 000 000000 000 000 000 000000 000000000...

----- 808 0000 000000 000000nm 1064 808 755 00nm 0 00000000 0000000 000000 00 0000.

---- 1 00000 000000000 .0000 10-1 0000000 .0000 000000 000000 !!! -

NEW



NEW



NEW



.میتواند به شما کمک کند تا به راحتی و بدون درد به نتیجه دلخواه خود برسید
میتواند 24 ساعته به شما کمک کند تا به نتیجه دلخواه خود برسید!

Free to be yourself

1 minute of painless hair removal,
Feeling the comfort and coolness that
has never been seen before.

14-inch screen

2 year warranty + work 24 hours



000000 000

350g lightweight **handle**, more ergonomic

Available in **black or white**





} 35mm



before after

Large 3.5cm spot size, and integration of three wavelengths
Hair removal is quick, efficient and convenient







One minute cooling -21°C

Best **TEC Cooling** system
hair remover laser



laser rod imported from America. 10 million—40 million shots. Permanent hair removal




755 808 1064nm 10000000 40000000 10000000 10000000 10000000.

* 100000000 100000000 100000000 100000000 100000000 100000000 100000000

(10000 10000 10000 10000) 10000 100000 40-10 *

10 *bar 12bar 16bar 1000000 10000000 100000000.

* 1600 / 1200 / 1000 1000000 1000000



1000000 1000000

Before and After

Treatment comparison



