



**MULTI-FUNCTION  
4 IN 1 MACHINE**

Produktbeschreibung

# Real shot display



808 handle



Elight handle



ND YAG Laser handle



RF handle



## **DiodenlaserElight/IPLNd yagRF, 4 verschiedene Griffe in einer Maschine**

1. Dieser 4 in 1 RFElight IPLDiodenlasernd Yag-Laser hat hochwertiges Zubehör, starke Leistung\*  
Diodenlasergriff verwenden Deutschlandimportierten goldenen Laserbarren, 50 Millionen Schüsse  
Zeit.


\*Elight-Griff verwendet UK-Lampe, 150.000-30.000 Aufnahmen

\* Nd yag: Rüsten Sie 5 Behandlungsköpfe aus, einstellbar 1064 nm/532 nm, fest 1064 nm/532 nm,  
fest 1320 nm.


\*RF: Rüsten Sie 6 Behandlungsköpfe aus, die zur Faltenentfernung, Hautstraffung und  
Hautstraffung verwendet werden.

## *Three wavelengths work together*

**12\*24mm<sup>2</sup>** {




**12\*35mm<sup>2</sup>** {



**Suitable for hair removal on any part of the body,  
Such as arms, legs, face, armpits, back, chest, etc.**





# ELIGHT HANDLE

## specification

- High-quality: Elight handpiece with super cooling system
- Optional Spot Size 15\*50mm, 15\*35mm
- Long shots time: 150000 times
- 7 filters: 420/510/560/585/640/690/755nm  
Standard: 420/510/680nm



## Before and after



Skin rejuvenation



Freckle removal



wrinkle reduction



Acne clearance



Hair removal



Pigmentation treatment

# HANDLE-ADVANTAGE



## specification

- Multiple choice  
fixed 1064nm & 532nm, adjustable 1064nm & 532nm (adjust 1-7 mm) for tattoo removal, nevus removal etc.
- Multifunction  
1320nm for skin rejuvenation
- Indicator light  
Dot position is the indicator light position dot more accurate



Adjustable  
532nm



Adjustable  
1064nm



Black doll  
1320nm



Fixed  
532nm



Fixed  
1064nm



# Before and after



## RF Handl e

12mm eye-rf	28mm small tripolar rf head	36mm big tripolar rf head	14mm bipolar-rf	28mm bipolar-rf	38mm bipolar-rf

# BEFORE & AFTER



Hair removal



Skin rejuvenation



Freckle removal



Acne removal



Vascular



Lose weight

## Advantage



### 808 handle

Painless hair removal, Permanent hair removal, Suitable for all skin tones.

### Elight handle

Elight for hair removal, skin rejuvenation, vascular therapy, pigmentation therapy, acne therapy.

### ND YAG Laser handle

ND Yag laser for tattoo removal, eyeliner removal, eyebrow embroidery removal, freckle removal.

### RF handle

Skin tightening/lifting, wrinkle removal on body

**Für weitere Informationen zur Produkteinführung und Verwendungsdetails wenden Sie sich bitte umgehend an uns Online 24 Stunden!!**