



¡13 años de fábrica profesional!

¡10 años de proveedor dorado de Alibaba!

¡5,0 alto punto de retroalimentación de Alibaba! ¡Servicio ODM/OEM!



14 inch
Screen Size



3000W
Input power



1-10HZ
Frequency



755nm 808nm 1064nm
Wavelength

>>Machine 14-inch large screen

>>3000W machine power to ensure the treatment effect

>>1-10Hz fast frequency and 15x25mm spot size ensures faster, smooth and effective treatment

>>808 single wave or 755nm/808nm/1064nm perfect triple wavelengths

>>Easy to operate Data Base system,can do the treatment by man mode and women mode, también diferentes partes del cuerpo. Para garantizar un parámetro seguro para el nuevo operador

>> Long life time More than 40 millions shots

Módulo de diodo láser dorado, Anti alta temperatura.



>>Double filters to keep water quality

>>TEC cooling Air Cooling Water CoolingSemiconductor cooling

, mantenga la temperatura del agua entre 26 °C y 28 °C todo el tiempo

>>24 hours continues work

>>Easy to operate Data Base system,can do the treatment by man mode and women mode,


también diferentes partes del cuerpo. Para garantizar una caja fuerte parámetro para el nuevo operador

>> Long life time More than 40 millions shots

Módulo de diodo láser dorado, Anti alta temperatura.



Descripción del Producto



El mango pesa solo 350 g, lo cual es conveniente para el tratamiento.

>>12 pieces or 16 pieces Germany Dilas original laser bars,600/800W/1200W power is capable to meet all demands of hair removal on the whole body;

>>1-10Hz fast frequency and 15x25mm spot size ensures faster, smooth and effective treatment;

>>Super sapphire touch cooling system ensures painless and comfortable treatment.

1) La pantalla puede agregar su propio logotipo

2) pantalla con 6 idiomas

3) diferentes partes del cuerpo para elegir, iun toque puede comenzar su tratamiento!

CERTIFICATE

