

# Large Spot Size

**Quick and Comfortable** 

With large spot size of 1.2\*3.5cm<sup>2</sup>, for faster treatments and better patient experience.







### **Continuous Contact Cooling**

TEC is an advanced technology that continuously cools skin and facilitates range-controlled temperature throughout the entire treatment. Patients can now enjoy a completely pain-free, cool and comfortable experience throughout the session.

### Manejar

- \*Barra láser de oro importado, Disparos hasta 40 millones
- \*Opciones de alta potencia de 10bar 1000W, 12 bar 1200W y 16 bar 1600W
- \*Tamaño de punto 12\*20mm, 15\*25mm, 12\*38mm para la opción
- \* Mango de 450 g, peso ligero, tratamiento de deslizamiento rápido gratuito
- \* Con luz indicadora, para mostrar las condiciones de trabajo (en la máquina y el mango), uso seguro


## **Handle Details**

### Three Wavelengths Work Togher




### Uso seguro en todo tipo de piel.

Alex 755nm para pieles claras Speed Diode 808nm para piel amarilla Long Pulse 1064nm para pieles más oscuras



- 755nm for white skin (fine, golden hair removal )
- 808nm for yellow neutral skin (brown hair removal)
- 1064nm for black skin (black hair removal)

Efecto del tratamiento



### Programa de la máquina

- 1. Pantalla grande de 12,4 pulgadas, toque rápido
- 2. Hombre y mujer, 6 partes del cuerpo, 6 tipos de piel que configuran software inteligente
- 3. 6 opciones de idioma, se puede agregar su idioma nativo
- 4. Servicio de programa OEM gratuito (logotipo/sistema...)