

Stronger, Faster

×

14-inch HD 4K Screen



超音波診断装置 -- V23

画面サイズ **14** インチ 画質 **HD 4K**

超音波診断装置 超音波診断装置 超音波診断装置

超音波診断装置 超音波診断装置 超音波診断装置

NEW



NEW





00000000 0000 00000 0000 000 00000000 000000 000000 000000 00 000000000 000000 0000 000000
0000...

1 0000 00000 0000 .0000 10-1 000000 .0000 00000 000 00000!!!-

350g lightweight **handle**, more ergonomic

Available in **black or white**



One minute cooling -21°C

Best **TEC Cooling** system
hair remover laser



1064 808 755 nm 1000 1200w/1600w 10 12bar 16bar 40-10 *

*. 1000 1200w/1600w 10 12bar 16bar 40-10 *
(1000 1200w/1600w 10 12bar 16bar 40-10 *)

* 1000 1200w/1600w 10 12bar 16bar 40-10 *

Before and After

Treatment comparison



.0000 0000 00 00 0000000 00000 00000 000000 00000 000000 0000
00000 00000 **24** 0000000 0000000 00000 00!