



# Diode Laser

## Permanent hair removal

[Understand more](#)

13 吋 大螢幕 操作簡便 易學!  
 10 年 經驗 專業 技術 團隊 服務 全球!  
 5.0 星 評價 ! 品質 保證 價格 公道 歡迎 ODM/OEM!

<p>14 inch Screen Size</p>	<p>3000W Input power</p>	<p>1-10HZ Frequency</p>	<p>755nm 808nm 1064nm Wavelength</p>

>>Machine 14-inch large screen  
 >>3000W machine power to ensure the treatment effect

>>1-10Hz fast frequency and 15x25mm spot size ensures faster, smooth and effective treatment

>>808 single wave or 755nm/808nm/1064nm perfect triple wavelengths

>>Easy to operate Data Base system,can do the treatment by man mode and women mode,

□□□□ □□□□□□ □□□□ □□□□ □□□□□□ □□ □□□□□□ □□□□□ □□□□ □□□□ □□□□□□

>> Long life time More than 40 millions shots

□□□□ □□□□ □□ □□□□□□ □□□□ □□□□□ □□□□□.



>>Double filters to keep water quality

>>TEC cooling Air Cooling Water CoolingSemiconductor cooling

□ □□□□□ □□□ °C 28-26 □□□□□ □□ □□ □□□□

24<< hours continues work

>>Easy to operate Data Base system,can do the treatment by man mode and women mode,

□□□□□ □□ □□□□□□□ □□□□□ □□□□ □□□□ □□□□□□□□

□□□□ □□□□□□□□ □□□□ □□□□□□□□


>> Long life time More than 40 millions shots

□□□□ □□□□ □□ □□□□□□ □□□□ □□□□□ □□□□□.



- Safety
- Quick
- Cooling system
- Permanent hair removal
- Best thermostatic water cycle system

□□□□ □□□□□□



□□□ □□□□ □□□□ □□□□ □□ □□□□ □□□ □□□ □□□□ □□□□ □□□□ 350 □□□□ □□□□

>>12 pieces or 16 pieces Germany Dilas original laser bars,600/800W/1200W power is capable to meet all demands of hair removal on the whole body;

>>1-10Hz fast frequency and 15x25mm spot size ensures faster, smooth and effective treatment;

>>Super sapphire touch cooling system ensures painless and comfortable treatment.

1) □□□ □□□□□ □□ □□□ □□□ □□□□□ □□ □□□□

□□□□ 6 □□ □□□□□ □□□□ (2

□□□ □□□□ □□ □□□ □□□□□ □□□□□ □□ □□□ □□ □□□□□□□□ □□□□ □□□ □□□□□ □□□ □□□□□ (3!

# CE and FDA from TUV

## CERTIFICATE

**U.S. FOOD & DRUG ADMINISTRATION**

510(k) Premarket Notification

Device Classification Name: **External Laser Unit**  
 510(k) Number: **K20025**  
 Device Name: **Medical Grade Laser**  
 Applicant: **Weifang Mingliang E No 13A26, Yuyang 9 Weifang, CN 26104**  
 Applicant Contact: **Li Mingliang**  
 Correspondent: **Beijing Believe-Med Rm 912, Building F, Fengshan District Beijing, CN 102401**  
 Correspondent Contact: **Ray Wang**  
 Regulatory Number: **EC 833**  
 Classification Product Code: **GE2**  
 Date Received: **03/02/2020**  
 Decision Date: **05/27/2020**  
 Decision: **Substantially Equivalent**  
 Regulation Medical Speciality: **General & Plastic S**  
 510k Review Panel: **General & Plastic S**  
 Type: **Traditional**  
 Reviewed By Third Party: **No**  
 Combination Product: **No**

**Certificate of Compliance**

No. 0417003-0419C02  
 Technical Construction No. 04190170185/0402

Certificate's Holder: **Weifang Mingliang Electronics Co., Ltd**  
 Ray 3, Workshop 2, No.13A26, Yuyang 8th Street, Yuyang Community, Shichang Street, High-tech Zone, Weifang, Shandong, China

Certification ECM Mark:

Product: **LASER BEAUTY DEVICE**  
 Model(s): **(See the following annex)**

Verification to: **Standard: EN 60601-2-22:2012 EN 60601-1:2006+A2:2014**  
 related to CE Directives: **Medical Devices Directive 93/42/EEC & 2007/45/EC**

Remarks: The products have been verified in a written report. The certificate certifies the achievement of the Certificate Mark of CE in accordance to the above listed standards. The above Certificate mark can be affixed on the products according to the CE regulation about its model and its use. The application procedure is found at [www.ecmcert.com](http://www.ecmcert.com). The Certificate of Compliance can be checked for validity in our website [www.ecmcert.com](http://www.ecmcert.com) upon payment of the production of the products.

Additional Information: **Application issued to CE mark**  
 The certificate is valid for 10 years. The CE marking procedure is in 2016. However, the certificate is a confirmation to use the CE marking Certificate Mark through an approved holder. Please use the product of the holding activities, and ensure that the structure and components be maintained, before using the CE mark on the products.

Date of Issue: **22 December 2019** Expiry Date: **22 December 2029**

CE Marking Authority:

ICM  
 11, Welterstrasse 5 33104 Bielefeld, Germany  
 www.tuv.com www.tuv.cn

**ITC**  
 Management Systems Certification Body No. 3022  
 Institute for Testing and Certification, Inc.  
 1, Bell 206, 754 27 Zia, Czech Republic  
 www.itc.com

**CERTIFICATE**  
 No. 10 0144 SJ

We certify on the basis of a performed audit that company  
**Weifang Mingliang Electronics Co., Ltd.**  
 13# Weifang Feng Cheng Street Weifang China  
 Company Reg. No.: 3707610000048

has implemented and documented a functional quality management system  
 in compliance with the requirements of the standard

**EN ISO 13485:2003**

Covering the following activities:  
**Design & Development, Production and Distribution of IPL, Phototherapy System, MD, YAG Q Switched Laser and RF Device**

The Certificate is based on the basis of the results mentioned in Audit Report No. 2019-08-00213. The Certificate shall be continued by positive results of surveillance audits, which the certified company continuously maintains.

Date of Issue: **05.01.2020**  
 Valid until: **05.01.2024**

By:

1st. Person Name: **Paul Vlas**  
 Head of Certificate Body

