

# Diode Laser

ice Hair removal



1800 1200 10-1

1. .
  2. : 36 \* 12 24 \* 12
  3. : 10-1
- 1000000  
1000000 1000000  
1000000 1000000



Medical 1200W/1800W

50 Million Shots

USA / Germany Imported Laser Bar



755nm+808nm+1064nm



ALEX 755nm

**For white skin**

(fine, golden hair removal )



SPEED 810nm

**For yellow neutral skin**

(brown hair removal)



YAG 1064nm

**For black skin**

(black hair removal)


□□□ 1800 □□□ 15 □□□ 1200 □□□ 10 □□□□ □□□□:

\* (□□□□□ □□□□ □□□ □□□□ □□ □ □□□□□□□□□□) .□□□□□ □□□□□□ □□□□□□ □□□□□

□□□ □□□□□□□ 50-20 □□□□ :□□□□ \*

1064 808 755 \*nm □□□ □□ □ □□□□ □□□□ □□□□□ □□□□□ □□□□□





00000 0000 00000

1) .000 00 00000000 00000 14 0000 0000 000000 0000 :0000 00000000

.0000 0000000 000000 000000 0000 00 :00000000 0000 (2

000000 (3 OEM: 00 0000 0000 00000 00000 0000000 00000 0000 0000000 00000 00 00000 00000 00 00000 00 00 0000 0000 00000000.