



.00000 0000 00000 0000 000000
00000 000 00 0000 000000 00000000 00000 *

00000 000 00 0000 000000 *

00000 0 00000 *

.00 0000 00000000 0000 000000 000000000000 *

808 000000 00000 0000 0000 *nm 7558081064 00nm.

350g

Ultra light handle



12×36mm

Handle

1800W diode laser hair removal




12×24 mm

Handle

1200W diode laser hair removal

이 레이저는 피부의 멜라닌 색소를 표적으로 하여 1064nm, 808nm, 755nm 파장의 레이저를 사용하여 피부의 멜라닌 색소를 파괴하여 모발을 제거합니다.

755nm (Alexandrite 레이저)는 피부의 멜라닌 색소를 표적으로 하여 모발을 제거합니다.
808nm (Diode 레이저)는 피부의 멜라닌 색소를 표적으로 하여 모발을 제거합니다.
1064nm (YAG 레이저)는 피부의 멜라닌 색소를 표적으로 하여 모발을 제거합니다.






ALEX 755nm
For white skin
 (fine, golden hair removal)




SPEED 810nm
For yellow neutral skin
 (brown hair removal)



YAG 1064nm
For black skin
 (black hair removal)



Treatment Comparison





Emergency
stop button

12 inch big screen
quick touch
easy operate

Diode laser handle

Vents

Plug

12 inch big screen quick touch easy operate