

1. .00 0 0000 00000 0000 00000 0000 00000 0000
 000 0000 36 * 12 0000 0000 0000 24 * 12 0000 000000 :0000 .2
 000 000 0000 00000 000 0000 00000 00000 10-1 000000 :0000 00000 .3.
 00000000 00000000





上海圣兰西医疗美容有限公司
上海圣兰西医疗美容有限公司



.....

..... 24*1.12
36*2.12mm

Arm hair removal



Face hair removal



Lip hair removal



Leg hair removal



Back hair removal



Underarm hair removal



Bikini hair removal





1064nm 808nm 755nm : 1064nm 808nm 755nm

755nm . 808nm . 1064nm .
 808nm . 1064nm .
 1064nm .

Medical 1200W/1800W

50 Million Shots

USA / Germany Imported Laser Bar



755nm+808nm+1064nm



ALEX 755nm

For white skin

(fine, golden hair removal)



SPEED 810nm

For yellow neutral skin

(brown hair removal)



YAG 1064nm

For black skin


(black hair removal)

1800 15 1200 10 :

* () .

50-20 :

1064 808 755 *nm



000000 00000 000000

1). 0000 00 00000000 000000 14 00000 00000 000000 00000 :000000 00000000

.00000 00000000 0000000 000000 0000 00 00 :000000000 00000 (2

000000 (3 OEM: 00 0000 0000 000000 000000 00000000 000000 0000 00000000 00000 00 00000 00 00 00000 00 00 0000 0000 00000000.

