



Medical Diode Laser Hair Removal

1800W+1200W handles

1. 100% 100% 100% 100% 100% 100% 100% 100% 100% 100%

1. 100% 100% 100% 100% 100% 100% 100% 100% 100% 100%
2. 100% 100% 100% 100% 100% 100% 100% 100% 100% 100%
3. 100% 100% 100% 100% 100% 100% 100% 100% 100% 100%
4. 100% 100% 100% 100% 100% 100% 100% 100% 100% 100%
5. 100% 100% 100% 100% 100% 100% 100% 100% 100% 100%

100% 100%

12*24mm or 12*36mm handle

Suitable for any unwanted hair on areas like face, arms, armpits, chest, back, bikini, legs...



- >>> Big 20*12 □□□□ □□□□□mm, 12*35mm for option
- >>> Laser technology, big power 10bar 1200W, 12 bar 1800W
- >>> 300g handle, light weight, free quick slide treatment



>>> New nose hair treatment head, applicable to small part treatment
>>> 1 minute refrigeration frosting, painless treatment experience






ALEX 755nm
For white skin
(fine, golden hair removal)



SPEED 810nm
For yellow neutral skin
(brown hair removal)




YAG 1064nm
For black skin
(black hair removal)



Treatment Effect

This is the clinical result of hair removal
It takes 3--5 times to realize
permanently hair removal





- * 4 ساله گارانتی و 14 ساله خدمات پس از فروش رایگان.
- * سیستم OEM برای بیمارستانها و کلینیکها.
- * قیمت بسیار مناسب و قابل پرداخت به صورت اقسامی.

مشخصات از دستگاه لیزر دایود موهای زائد

مورد	755nm/808nm/1064nm لیزر دایود عمودی
نوع لیزری	نوار لیزر پزشکی لیزر دایود: 10 بار / 15 بار
قدرت هر نوار	120w/بار پزشکی
سیستم نمایش	صفحه نمایش لمسی 14 اینچی رنگی واقعی
طول موج	808nm/755nm/1064nm
Spot Size	24*12mm ² /12*36mm ² پزشکی
طول عمر دسته لیزری	حداقل 40 میلیون شات
انرژی	(اختیاری cm ² /قابل تنظیم) تا 150 J / cm ² 10-150
مدت زمان نبض	میلی ثانیه قابل تنظیم 30 ~ 400
فرکانس	هرتز؛ 1.0-10
سیستم خنک کننده	کمپرسور هوا گردش آب را ببندید
قدرت دسته	1000w/1200w/1600w/1800w
قدرت ماشین	4000 وات
منبع برق	ولت / 50 هرتز یا 110 ولت / 60 هرتز 220
دمای کریستال	تا 3 درجه سانتی گراد -20
دمای ماشین	درجه سانتی گراد 24
نیروی مورد نیاز	AC220V/110V 50~60Hz

□□□□ □□□□

Free Shipping

Custom Logo

Hot Selling Big Spot



12×36mm

Medical diode laser Hair removal machine

V19



• 100% 100% 100% 100% 100% 100%
100% 100% 100% 100% 100%
100% 100% 100% 100% 100% 100% 100% 100% 100% 100%
100% 100% 100% 100% 100%
100% 100% 100% 100% 100% 100% 100%
• 100% 100% 100% 100% 100% 100% 100% 100%
100% 100% 24!