



000 0000 000000 000000

0000 0000 5.5 = 00000 30 0000 00000 0000 .1

000 0000 00000000 0000 000000 000 00000 .2

0000 00000 19 0000 .3

000000 0000 00000 16 000000 .4

0000 0000 .5

0000 00 000 00 000000 000000 .6

000000 0 0000 0000 000000 .7

000000 000000 0000 0 0000 0000000 .8



Operation



Hamstring



Thighs



Abdomen




Arm



Hip



BEFORE AND AFTER



Before

5 Weeks After

8 Weeks After



Before

After



Before

After







15 Years Laser Beauty Device Manufacturer

Make your beauty business easier



Smayer is a professional manufacturer of laser beauty devices. We have 15 years of experience in the industry. Our products are widely used in beauty salons and clinics. We provide OEM and ODM services for our clients. We have a complete production line and strict quality control system. We are committed to providing high-quality products and services to our customers.

Our Factory Show

**15 Years
Professional Factory**




**Professional
Production R & D
Sales Department**



**Mass production
of various styles**



**Real Feedback
of Alibaba**



Certificate display

Certificate is the best guarantee for you



□□□ □ □□□ □ □□□□ □□□□

Packing

Safe and thoughtful packaging service



Tight packaging

The machine shall be made of waterproof materials and packaged

Anti collision foam board inside aluminum box



3-5 business days delivery

Professional guide:

24 Hours Professional online customer service

Our company has a professional packing team, 7-3:00 AM-11:00 PM, we can provide you with professional packing service.
We can provide you with professional packing service, DHL, TNT, FedEx, UPS, etc. We can provide you with professional packing service.
We can provide you with professional packing service, 24 hours online customer service, we can provide you with professional packing service.

CONTACT US



In terms of following trade details:

- 1)The shipping cost are based on DHL cost to your country
- 2)Delivery time:3-5 days after receiving payments.
- 3)Payment terms:T/T(bank transfer),Western Union.
- 4)Warranty:2 years for the whole machine
- 5)Warranty Policy:one month guarantee replacement (we bear all shipping cost),2years warranty,lifelong maintenance.
- 6)Packing method: Aluminum box with hard plastic foam.
- 7)Free-training to customers:Free We could offer DVD and User's manual and 10 Hours online service to our customers
- 8)Certifications:CE TUV SFDA ISO 13485
- 9)MOQ:1 UNIT
- 10)The quotation here above is valid within 1 month

問問問問 問問問問

Q1. 問問 問 問問 問問問 問 問 問問問問 問

A: 問問 問問問問 問問問問 問問 7 問問 問 問問問問問 問 問.

Q2. 問問問 問 問問問 問問問問 問問問 問問

A: 問問 問 問問問 問問 問問 問 問問問問問問 問問問問 問 問 問問問 問 問問問 問 問問問問 問問問 問問問 問.

Q3. 問問 問問 問問 問問問問 問問

A: 問問 問問問 問問問 問問問問 問問問 問 問問問 問問 問問 問 問問問 問問問 問 問問問問問 問問.

問問 問問問 問問 問問問 問 問問 問.

Q4. 問問問 問問問 問 問問問問 問問問

A: 問問問 問問問 問問問問 問 Alibaba 問 問問問問 Western Union VISA 問問問問 問 問.

Q5. 問問問 問問問問

A: 2 問問 問問問問 問 問 問問問問問 問問問 問問 問 問問問 問問 .問問 問問 問問 1 問問問問 問問 問問 問問 問問 問 問問問 問問 問問.