



1064/808/755 nm laser treatment for skin rejuvenation. 1 or 3 sessions.

* Laser treatment for skin rejuvenation. 1 or 3 sessions. 808 nm laser treatment for skin rejuvenation. 1 or 3 sessions. 755 nm laser treatment for skin rejuvenation. 1 or 3 sessions.

Suitable for hair removal on any part of the body,
Such as arms, legs, face, armpits, back, chest, etc.



• Cooling effect •

The handle is frosted for one minute,
making your treatment more comfortable and painless






Soon

Permanent hair removal

1064nm YAG, 808nm Speed, 755nm Alex, 4mm treatment head.

755nm Alex (epidermal/dermal treatment) for epidermal and dermal hair removal.
 808nm Speed (epidermal) for epidermal hair removal.
 1064nm YAG (dermal) for dermal hair removal (1064nm).






ALEX 755nm
For white skin
(fine, golden hair removal)




SPEED 810nm
For yellow neutral skin
(brown hair removal)



YAG 1064nm
For black skin
(black hair removal)



Treatment Comparison





Emergency stop button

12 inch big screen quick touch easy operate

Diode laser handle

Vents

Plug

1. 2. 3. 4. 5. 6.

1. .
2. .
3. .
4. .
5. .
6. .

!