

# Medical



1200W+1800W high power handle  
14-inch 4K HD large screen  
755nm+808nm+1064nm  
Painless permanent hair removal

# Diode Laser

□□□□ □□□□ □□□□ □□□□ □□□□ □□□□ □□□□ □□□□

4 □□□□ 14 □□□□ □□□□ □□□□K HD

1200 □□□□ □□□W1800W □□□□ □□□□ □□

□□□ □□□□□□□□□□ □□□□ □□ □□□□□□□□ □□□□□

□□□□ □□ □□□ □□□□□□□□ □□ □□□□ □□□□ 350 □□□□ □□□□ □□□□

□□ □□□□□ □□□ □□□□□ □□□□□□□□ □□□□ □□□ □□ □□□□ □□□□ 36 \* 12 □□□□ □□□□ 24 \* 12 □□□□□

□□□ □□ □□□□□□□□

# Machine Display



1. .0000 00000 000 0000 000 0000 000 00 00000 0 00 0000000000 000 0000 000 000 00000000 0000000000 1064 0000000000 808 0000000000 755 00000 000000 00000 00000000 .
- 2.
- 3 000000 000 000000.TEC 000000 1 00000 000000 0000



# Before and After

## Treatment comparison



□□□□ □□ □□ □□□□ □□□□□□ □□□□ □□□□ .□□ □□□ □□□□ □□□ □□□□ □□□□ □100 □□□□□□ □□□□!