

7D Focused Ultrasound



SHURINK L4-4.5
SHURINK DL4-5

FACIAL SCULPTURE

Smaller, more convenient, more accurate, more suitable for facial skin

- 1.5 mm
- 3.0 mm
- 4.5 mm

7D HF HANDLE

More comprehensive, more detailed and more suitable for comprehensive treatment of the body


- MF2
- MF6
- MF9
- MF13
- L4-4.5
- L7-1.5
- L7-3.0



EXFU is a professional medical device for facial and body treatment. It is designed to provide a comprehensive, detailed, and accurate treatment for facial skin and the body. The device features a large screen for monitoring and control, and multiple handles for different treatment areas. The handles are labeled with their respective treatment parameters, such as 1.5 mm, 3.0 mm, 4.5 mm, MF2, MF6, MF9, MF13, L4-4.5, L7-1.5, and L7-3.0. The device is suitable for a wide range of skin types and conditions, and is used by professional medical staff. The device is also suitable for home use, and is a valuable addition to any medical or beauty practice. The device is designed to be safe, effective, and easy to use, and is a reliable choice for medical and beauty professionals.




EXFU




Face Treatment



Body Treatment



Body Treatment



• Treatment Effect •



000	000000 0000 7 000000 00 000000 0000000000
000	0000 000 00 000000 0000000000

探头 探头频率	4.5 探头 3.0 探头 1.5 :探头
探头 探头频率	7D HF: MF2/6/9/13 L4-4.5 L7-1.5 L7-3.0
探头 探头频率	0.1J-2.5J (0.1J 探头 探头)
探头 探头 探头	1.0mm-4.0mm
探头 探头	5-25 探头 探头
探头 探头	探头 300000 探头 探头
探头 探头 探头	15 探头
探头	110V-240V 50/60Hz



探头 探头 探头 探头 探头 探头 探头
 探头 探头 探头 15.6 探头 探头 探头 探头 .1
 .探头 探头 探头 探头 6 探头 探头 探头/探头 .2
 .探头 探头 探头 探头 .3
 探头 探头 .4 OEM 探头/探头) 探头...)
 探头 探头 探头 探头 探头 探头 探头!