



Medical Diode Laser Hair Removal


1800W+1200W handles



00000 0000000



0000 00000






Laser show Handle laser bar optional



④ 10bar

④ 12bar

④ 16bar



□□□□ □□□□□□ □□ □□ □□□□

Suits for all skin and hair types



808nm

1064nm



- **755nm for white skin** (fine, golden hair removal)

● **808nm for yellow neutral skin** (brown hair removal)

- **1064nm for black skin** (black hair removal)




- 1.** _____ **15** _____.

2-□□ □□ □□□ □□□ □□ □ □□□ □□ □□ □□□□ □□□□ □□□ □□□□.

[illegible]

Treatment Comparison



□□□□ OEM

1. □□□□□ □□□ □ □□□
2. □□□□□□□ □□□□ □□□□□ □□ □□□□□ □□□ □□□□

customer visit



Packing

Safe and thoughtful packaging service



Tight packaging

The machine shall be made of waterproof materials and packaged

**Anti collision foam board
inside aluminum box**



3-5 business days delivery

Professional guide:

24 Hours Professional online customer service