



Macro & Micro High Intensity Focused Ultrasound System

7 CARTRIDGES + 4 GUN HEADS


1.5mm, 2.0mm, 3.0mm, 4.5mm for face treatment
6.0mm, 9mm, 13mm for body treatment.




MF2.0MM

Unique ultra-thin probe Clexis patented technology


"Special cartridge for eyes treatment"



MF2.0MM is a non-invasive skin treatment device that uses a unique ultra-thin probe to deliver energy to the SMAS layer. The device is designed for eye treatment and is safe and effective. The energy is delivered in a controlled manner, ensuring that the treatment is precise and targeted. The probe is made of a special material that allows it to penetrate the skin without causing any damage. The energy is converted into heat, which stimulates the production of collagen and elastin, resulting in tighter and smoother skin. The treatment is suitable for all skin types and ages. The device is easy to use and can be performed by a trained professional. The results are visible immediately and last for several months. MF2.0MM is a revolutionary skin treatment device that offers a safe and effective solution for eye treatment.



Body Treatment



• *Treatment Effect* •



□□□	□□□□□□ □□□□ 7 □□□□□□ □□ □□□□□□ □□□□□□□□□□
□□□	□□□□□ □□□ □□ □□□□□□ □□□□□□□□□□

□□□□ □□□□□□	□□□□ □□□□ 4.5 □□□□ □□□□ 3.0 □□□□ □□□□ 1.5 :□□□□
	7D HF: MF2/6/9/13 L4-4.5 L7-1.5 L7-3.0
□□□□□□ □□□□□□	0.1J-2.5J (0.1J □□□□□□ □□ □□)
□□□□ □□□ □□□□□□	1.0mm-4.0mm
□□□□ □□□□	5-25 □□□□ □□□□□□
□□ □□□□	□□□□ 300000 □□□□□□□□ □□
□□□□□□ □□□□□ □□□□□□□□	15 □□□□□□
□□□□□□	110V-240V □ 50/60Hz

Comparison before and after use



- □□□□ □□□□□ □□□□ □□□□ □□□□□□□□
- □□□□□ □□□□□□ 15.6 □□□□ □□□□ □□□□□ □□□□ .1
- .□□□□□□ □□□□ □□□□ □□□□ 6 □□□□ □□□□□ □□□□/□□□□ .2
- .□□□□ □□□□□ □□□□ □□□□□□□□ .3
- □□□□□ .4 OEM □□□□□□/□□□□□□) □□□□□□□□□□...

EXFU Seiko



11D HIFU

Macro & Micro High Intensity Focused Ultrasound System



SAFE



EFFECTIVE



MULTI-DEPTH



□□□□ □□ □□□□ □□□□ □□□□□ □ □□□□ □□□□!