



በቅድመ ማግኘት ለሚገቡ ጉዳዮች ማግኘት፡

- 1) በቅድመ ማግኘት ማግኘት ማግኘት ማግኘት.
- 2) በቅድመ ማግኘት ማግኘት 168 ሰዓት ማግኘት ማግኘት
- 3) በቅድመ ማግኘት ማግኘት ማግኘት ማግኘት ማግኘት ማግኘት ማግኘት
- 4) በቅድመ ማግኘት ማግኘት ማግኘት ማግኘት

# Machine Parts



|      |                                |
|------|--------------------------------|
| 产品名称 | 808 光子脱毛仪 A1                   |
| 品牌   | 808 光子脱毛仪                      |
| 规格   | 755 光子脱毛仪 1064 光子脱毛仪 808 光子脱毛仪 |
| 适用部位 | 4 光子脱毛仪 13.3 英寸 4K 真彩触摸屏       |
| 功率   | 1-300 瓦                        |
| 重量   | 1-10 公斤                        |
| 尺寸   | 14*14 厘米                       |
| 能量   | 1-150J/cm2                     |
| 电压   | 220V/110V 50HZ/60HZ            |
| 频率   | AC220V/110V 50HZ/60HZ          |
| 波长   | 355*390*422mm                  |
| 材质   | 铝合金                            |
| 包装   | 2 个                            |



# Machine Display



Left 45°



Positive



Right 45°

## Handle Display

*The handle is small and light, weighing as low as 168.8g.  
Very light to operate*



この装置は、治療時に患部を冷却する効果があり、痛みを軽減します。

24時間いつでも使用可能。冷却効果は持続的。冷却効果は持続的！

## Cooling Effect

*The handle has a cooling effect that improve treatment comfort.*



**Just One minute!**  
Quick cooling down to

**-25°C**

この装置は、治療時に患部を冷却する効果があり、痛みを軽減します。 :お問い合わせ 808

この装置は、治療時に患部を冷却する効果があり、痛みを軽減します。 :お問い合わせ 755

**ND YAG 1064** 装置は、治療時に患部を冷却する効果があり、痛みを軽減します。 :お問い合わせ



# Working Principle

755nm and 808nm and 1064nm , suits for all skin and hair types.



**ALEX 755nm**

**For white skin**

(fine, golden hair removal )



**SPEED 810nm**

**For yellow neutral skin**

(brown hair removal)



**YAG 1064nm**

**For black skin**

(black hair removal)

## Treatment Effect



Before



After



Before



After




Before



After



Before



After

**CERTIFICAT E**

BV-France's most authoritative  
certification body BV Group certification



TUV-Germany's most authoritative  
certification body Rheinland certification



# SGS

SGS-Switzerland's most authoritative  
certification body SGS certification



Italy ECM Certification



ISO  
13485

NQA medical system  
certification



!~~~~~ 0000 00 00 0000 00000000 0000 00 000000 00 0000 000 000  
000000 **24** 000000 000000!