



1. Three wavelength diode laser.
2. Golden laser bars from German or USA improved 40,000,000+ longer lifespan.
3. Safe, Fast, Painless, permanent
4. For all skin and hair style hair removal
5. 14.0 inch true color touch screen with 6 skin colors.
6. Two models can be selected, single handle or double handle



Most advanced TEC cooling system inside this machine to deliver a super good cooling results. 24 hours working time.Sapphire material.Actual Pain-free experience!

One minute cooling -21°C

Best **TEC Cooling** system hair remover laser



- 1)10 bars/ 12 bars/ 16 bars, 10-40 million shots.
- 2)Big spot 12*20mm / 12*35mm for fast hair removal.
- 3)Frozen ice cooling sapphire crystal, for painless hair removal.
- 4)Hight quality gold laser bar imporved from German.
- 5)Medical handle as low as 350g easy to treatment.

1 minute
and you're done

35mm {

FASTER

Screen system



1. 14 inch big touch screen, humanized designed.
2. Male/female, body parts, 6 skin types for option.
3. 6 language systems, native language + whatever you want.
4. Free OEM program service (logo/system...)

| | |
|-------------------------|--|
| Laser bar | Germany/USA imported, 8 bar, 10 bar, 12 bars, 16 bars |
| Wavelength | 755nm 808nm 1064nm/ or single 808nm |
| Energy | 0-150J/cm ² |
| Frequency | 1-10hz |
| Pulse width | 40-800ms |
| Operative LCD Interface | 14.0 TFT True Color LCD touch screen |
| Cooling system | Air and Close Water Circulation Treatment head with semiconductor cooling temperature adjustable |
| Power | 4000W |
| Electrical source | 220V/50HZ or 110V/60HZ |
| N.W/ G.W | 62kg/78kg |

**Standardized
production**



**Professional
Production R & D
Sales Department**



**Mass production
of various styles**



- >>OEM your own machine body.
- >>OEM your own program software.
- >>OEM pretty machine colors.
- >>Print your LOGO on machine body and screen.

Custom Section


CASE design



LOGO design



SYSTEM design



Certificate Display

CERTIFICATE DISPLAY

