

DIODE LASER



ICELLO



MLKJ Best commercial laser hair removal machine.

According to the principle of selective photothermodynamics, the laser wavelength, energy and pulse width can be properly adjusted to make the laser pass through the skin surface to the root hair follicle. Light energy is absorbed and converted into heat energy that destroys hair follicle tissue, thereby rendering hair unable to regenerate, making MLKJ lasers' patented technology painless and permanent. Laser hair removal device is a non-invasive treatment method, after treatment will not have a certain damage to the surrounding skin, will not affect the normal life and work.

☐ FDA & Medical CE Certified Machine

☐ 3 in 1 - Three Laser Technologies in one machine

☐ Alexandrite, diode and Nd:YAG lasers - treats all skin types

☐ Easy to operate

☐ Online training available

☐ 50 Million shots per handpiece!

Handle display



1. Pain-free. The tip of handpiece cools down to -32°C , makes treatment comfortable and effective

2. Early treatment results. 2-3 treatment sessions can show permanent results

3. Virtually pain-free contacted with sapphire incomparable comfortable, and safety

4. Facial Tip Improved Reach. Easily treat areas that are normally hard to reach, including the ears, nostrils and glabella.



755/808/1064nm laser hair removal machine price.

>> 808 single wave or 755nm/808nm/1064nm perfect triple wavelengths

>> America Imported Emitter, strong energy

>> 10/12 pieces or 16 pieces Germany Dilas original laser bars, 1200W/1800W power is capable to meet all demands of hair removal on the whole body


Medical diode laser

50000000+shots
can treat 10000+clients

- *Not easy to be burned
- *Low maintenance cost
- *Stable quality



- >>1-10Hz fast frequency,faster, smooth and effective treatment
- >>Super sapphire touch cooling system ensures painless and comfortable treatment
- >> Italy imported stainless water pump, it can keep machine working long time without stop. Make sure the heat were removed very quickly
- >>Double filters to keep water quality



The real experience is visible

Through live test verification ,safe ,long,painless
and effective



14-inch 4K screen, smart touch experience

- * 6 skin styles to choose from
- * Easy to operate providing operation skills/teaching services
- * Free logo and language can be customized on the screen



You'll get

Two major service guarantees
Worry-free after-sales

01 Professional after-sales group

24 hours offers the solution
Solve problems quickly with your active cooperation

02 Free maintenance

Warranty period, Send your replacement
parts for free

MLKJ

Manufacturer of beauty instruments for 15 years, have certificates FDA, CE ,ISO13485.
We are committed to making your beauty business easier!

Establishment Registration & Device Listing

FDA Home Medical Devices Databases



New Search Back To Search Results

Establishment:

WEIFANG MINGLIANG ELECTRONICS CO.,LTD.
Floor 3, Building 4,Second Accelerator, Yuqing Eastern Street, High-Tech District
Weifang Shandong, CN 261000
Registration Number: 3007549439
FEI Number*: 3007549439
Status: Active
Date Of Registration Status: 2022

Owner/Operator:

Weifang Mingliang Electronics CO.,LTD.
Floor 3, Building 4,Second Accelerator, Yuqing Eastern Street, High-Tech District
Weifang, Shandong CN 261000
Owner/Operator Number: 10082104

Official Correspondent:

Pingping Lv
Floor 3, Building 4,Second Accelerator, Yuqing Eastern Street, High-Tech District
Weifang, Shandong CN 261000
Phone: 86-10-56335780

US Agent:

Wen Xu
4611 Maydew Mews Way Ellicott City MD
Ellicott, MD US 21042
Phone: 443 8968583 Ext
Email: Infor@Believe-Med.Com

* Firm Establishment Identifier (FEI) should be used for identification of entities within the imports message set

Page Last Updated: 05/09/2022

Note: If you need help accessing information in different file formats, see Instructions for Downloading Viewers and Players.

Language Assistance Available: Español | 繁體中文 | Tiếng Việt | 한국어 | Tagalog | Русский | العربية | Kreyòl Ayisyen | Français | Polski | Português | Italiano | Deutsch | 日本語 | فارسی | English



Authentication

FDA,NQA13485,CB certification
and more Complete certification
materials



60000+ Customers

We have professional after sale
team and we provide 7*24
services



2 Years Warranty

2 years warranty for each
machine, All questions will get
reponse within 24 hours

Welcome to consult! 7X24 hours online service!