



## ***Hair Removal + ND YAG***

*2in1 Multifunctional combination machine  
one machine and two effects*

**Multi-functional beauty instrument manufacturers**

**2 in 1 Diode laser hair removal +nd yag laser tattoo removal**

Products Description

## 2in1 Multifunctional combination machine



- ✓ Hair removal
- ✓ Toner whitening
- ✓ Mole removal




- ✓ Tattoo removal
- ✓ Freckle
- ✓ Birthmark removal



# Real shot display

Higher-definition screen, stronger cooling effect,  
more effective treatment



## **Diode laser+Elight/IPL+Nd yag+RF, 4 different handles into one machine**

1.This 4 in 1 RF+Elight IPL+diode laser+nd yag laser have high quality accessories , strong power\*Diode laser handle use Germanyimported golden laser bar ,50 million shots time.


\*Elight handle use UK lamp, 150,000-30,000shots

\*Nd yag :Equip 5 treatment heads, adjustable 1064nm/532nm, fixed 1064nm/532nm, fixed 1320nm.


\*RF :Equip 6 treatment heads, used for wrinkle removal, skin lifting, skin tightening.

## *Three wavelengths work together*

12\*24mm<sup>2</sup> {



12\*35mm<sup>2</sup> {



Suitable for hair removal on any part of the body,  
Such as arms, legs, face, armpits, back, chest, etc.







# HANDLE-ADVANTAGE



## specification

- Multiple choice  
fixed 1064nm & 532nm, adjustable 1064nm & 532nm (adjust 1-7 mm) for tattoo removal, nevus removal etc.
- Multifunction  
1320nm for skin rejuvenation
- Indicator light  
Dot position is the indicator light position dot more accurate



Adjustable  
532nm



Adjustable  
1064nm



Black doll  
1320nm



Fixed  
532nm



Fixed  
1064nm



## Before and after



## 14inch true touch screen

6 kinds skin type for choice. accept ODM OEM Service.

The new **14-inch** screen is clear easy to operate.



For more information about product introduction and use details,

**please consult us immediately**

**Online 24 hours!!**