

## Product Description

---

### Advantages



# Diode Laser



15 Bar **1800W** Handle

\*Increase the treatment head of nasal hair  
and treat small parts easily

1. 3 waves **1064nm + Diode 808nm + 755nm**, 3 in 1 all skin type can be used  
High power: NEW 1000W/ 1200W /1600W handle for options
2. Big spot size **12\*20mm, 12\*35mm** face & body hair removal solution!
3. Cooling system: Unique Dual- TEC cooling system, most effective to ensure machine work all day, 24hr no down time.
4. Long life: double life support guarantee

### Handle of diode laser hair removal machine:



**Different spot size** Handle color optional



12\*20mm



12\*35mm



12\*20mm



12\*35mm

### Safe Use On All Skin Type

Alex 755nm for light skin

Speed Diode 808nm for yellow skin

Long Pulse 1064nm for darker skin

# Treatment effect



## Software & Program

1. Male & femal, 6 body parts, 6 skin types setting smart software
2. 6 language options, your native language can be added
3. OEM program service available(logo/system...)



Product name	Laser Ice platinum hair removal machine
Laser bar	Germany/USA imported, 8 bar, 10 bar, 12 bars, 16 bars
Wavelength	755nm 808nm 1064nm/ or single 808nm
Energy	0-150J/cm2
Frequency	1-10hz
Pulse width	